


A.B.S. Curved Back Assembly

Recline, Sliding Back (VSR)

Tilt/Recline and Recline Only



A.B.S. Curved Back Assembly

Recline, Sliding Back (VSR)

Tilt/Recline and Recline Only

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		0047800-NLA	NLA - Block, Glide	1	4
2		0052906-NLA	NLA - Screw, Socket Button Head (1/4-20 x 3/4")	1	4
3-6		See Chart	#'s 3-6 A.B.S. Curved Back Assembly	1	1
	A,E	1125723	Kit, Hook Fastener w/ Adhesive Back	1	1
	A,B	0143900-NLA	NLA - Flap, Back Pan - Tarsys	1	1
	A,C	0144000-NLA	NLA - Flap, Back Pan - Tarsys	1	1
	A,D	0144100-NLA	NLA - Flap, Back Pan - Tarsys	1	1
NOTE: A - Not Shown B - 16"-17" Wide C - 18"-19" Wide D - 20"-24" Wide E - Hook Fastener Adhesive Kit is used to attach the Back Cover to the Back Pan and must be ordered when a new Back Pan or Back Assembly is requested.					

Material Back Width x Height	Overall Back Height	Item #3 Back Pan	Item #4 1" Back Foam	Item #5 Back Cover	Item #6 Back Asbly w/Sliding Hdwr
14" X 14"	21"-22"H	0128819	1073068	1073089-NLA	1077095
15" x 14"	21"-22"H	0128919-NLA	1073069	1073090-NLA	1077096-NLA
16" x 14"	21"-22"H	0129019-NLA	1073070	0124821-NLA	1077097
17" x 14"	21"-22"H	0129119	0310820-NLA	0124921-NLA	1077098
18" x 14"	21"-22"H	0129219-NLA	1073071-NLA	0125021-NLA	1077099-NLA
19" x 14"	21"-22"H	0129319-NLA	1073072-NLA	0125121-NLA	1077100-NLA
20" x 14"	21"-22"H	0129419-NLA	1073073	0125221-NLA	1077101-NLA
21" x 14"	21"-22"H	0129519	1073074	1073091-NLA	1077102
22" x 14"	21"-22"H	0129619-NLA	1073075	1073092-NLA	1077103-NLA
23" x 14"	21"-22"H	Special	Special	Special	Special
24" x 14"	21"-22"H	Special	Special	Special	Special
14" x 16"	23"-25"H	0048819-NLA	1073070	0010821-NLA	1077104-NLA
15" x 16"	23"-25"H	0048919-NLA	1073077-NLA	0011021-NLA	1077105-NLA
16" x 16"	23"-25"H	0049019-NLA	0009720	0011221-NLA	1072878
17" x 16"	23"-25"H	0049119	0009920-NLA	0011421-NLA	1072879
18" x 16"	23"-25"H	0049219-NLA	0010120-NLA	0011621-NLA	1072880-NLA
19" x 16"	23"-25"H	0049319	0010320-NLA	0011821-NLA	1072881
20" x 16"	23"-25"H	0049419-NLA	0009820-NLA	0012021-NLA	1072882-NLA

A.B.S. Curved Back Assembly

Recline, Sliding Back (VSR)

Tilt/Recline and Recline Only

Material Back Width x Height	Overall Back Height	Item #3 Back Pan	Item #4 1" Back Foam	Item #5 Back Cover	Item #6 Back Asbly w/Sliding Hdwr
21" x 16"	23"-25"H	0049519-NLA	1073078	0304021-NLA	1077106-NLA
22" x 16"	23"-25"H	0049619	1073079	0164421-NLA	1077107
23" x 16"	23"-25"H	Special	Special	Special	Special
24" x 16"	23"-25"H	Special	Special	Special	Special
NOTE: Some foam size part numbers can accommodate two width and two back heights by turning the foam around 90 degrees. Example would be a (16" W x 20" H) can be (20" W x 16" H) in other applications.					