


Recessed Back Assembly

Sliding, VSR

Tilt/Recline & Recline Only



Recessed Back Assembly

Sliding, VSR

Tilt/Recline & Recline Only

Material Back Width x Height	Overall Back Height	Item #1 Back Pan	Item #2 Back Assembly w/ Sliding Hardware
14" x 15"	23"	0048019	0323800-NLA
15" x 15"	23"	0048119-NLA	0323900-NLA
16" x 15"	23"	0048219	0324000-NLA
17" x 15"	23"	0048319	0324100-NLA
18" x 15"	23"	0048419	0324200-NLA
19" x 15"	23"	0048519	0324300-NLA
20" x 15"	23"	0048619	0324400-NLA
21" x 15"	23"	Special	Special
22" x 15"	23"	Special	Special
23" x 15"	23"	Special	Special
24" x 15"	23"	Special	Special

Recessed Back Assembly Sliding, VSR Tilt/Recline & Recline Only

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-2		See Chart	#'s 1-2 Recessed Back Assembly (Sliding, VSR)	1	1
3		0047800-NLA	NLA - Block, Glide	1	4
4		0223506-NLA	NLA - Screw, Socket, Cap, Button Head (1/4-20 x 1/2")	1	4
5		0047902	Hook, Mounting	1	4
6		7930323-NLA	NLA - Package, Locknut (1/4-20)	25	8
7		0052906-NLA	NLA - Screw, Socket Button Head (1/4-20 x 3/4")	1	8
	A,B	0143900-NLA	NLA - Flap, Back Pan - Tarsys	1	1
	A,C	0144000-NLA	NLA - Flap, Back Pan - Tarsys	1	1
NOTE: A - Not Shown B - 16"-19" Wide C - 20"-24" Wide					