



AT5544 Legrest (Power)

Tilt Only or Tilt/Recline and Recline Only



*AT5544 rigging shown for reference only.
See 2G Storm Series™ parts catalog form
number 97-188 for legrest components.*

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		0025306-NLA	NLA - Push Pin, AT5544	1	1
2		0025602-NLA	NLA - Bracket, Push - Right	1	1
2		0025502-NLA	NLA - Bracket, Push - Left	1	1
3		0055206-NLA	NLA - Screw, Cap, Button Head, Socket (5/16-18 x 1-3/4")	1	2
4		0055900-NLA	NLA - U-Clamp (7/8")	1	1
5		0054806-NLA	NLA - Washer, Flat (11/32 x 11/16 x 1/16")	1	5
6		0054906	Locknut, High Profile (5/16-18)	1	2
7		0054106-NLA	NLA - Locknut, High Profile (3/8-16)	1	1
8	A	1077516-NLA	NLA - Hardware, Legrest Mounting, AT5554 - Right	1	1
8	A	1077517-NLA	NLA - Hardware, Legrest Mounting, AT5554 - Left	1	1

NOTE: A - Includes items 1-7