





# Pivot/Slide Tube and Footplates

## PH904A and PHAL4A

Heel loops are not included. See Footrest and Legrest Accessories.  
For complete pivot/slide tube and footplate assemblies, see Front Riggings Reference chart.

|   |   |   |
|---|---|---|
|  |  |  |
| 1350<br>Extra large   | 1450<br>Deluxe  | 1650<br>Standard Aluminum   |



# Pivot/Slide Tube and Footplates

## PH904A and PHAL4A

Heel loops are not included. See Footrest and Legrest Accessories.

For complete pivot/slide tube and footplate assemblies, see Front Riggings Reference chart.

| Item Number | Note  | Part Number | Description   | Quantity Per |          |
|-------------|-------|-------------|---|--------------|----------|
|             |       |             |   | Package      | Assembly |
| 1           |       | 1027097     | Plug Button, Black (3/4")   | 1            | 2        |
| 2           | A,B,D | 1100568-NLA | NLA - Pivot/Slide Tube w/ Plug Button, Black  | 1            | 1        |
| 2           | A,B,D | 1100569-NLA | NLA - Pivot/Slide Tube Assembly, Extra Large, Standard Aluminium - 2" Longer, Black | 1            | 1        |
| 2           | A,B,D | 1100570-NLA | NLA - Pivot/Slide Tube w/ Plug Button, Black (4" Longer)                            | 1            | 1        |
| 2           | A,B,E | 1041526     | Pivot/Slide Tube w/ Plug Button Assembly, Chrome                                    | 1            | 1        |
| 2           | C,D   | 1100571-NLA | NLA - Pivot/Slide Tube w/ Plug Button, Deluxe, Black                                | 1            | 1        |
| 2           | C,E,H | 1023804-NLA | NLA - Pivot/Slide Tube w/ Plug Button, Deluxe, Chrome                               | 1            | 1        |
| 3           | A     | 44200X027   | Footplate, Aluminum, Medium, Black (7-3/4" W x 6" L)                                | 1            | 1        |
| 3           | C     | 1043659-NLA | NLA - Footplate, Aluminum, Deluxe, Black - Right (8" W x 5-3/4" L)                  | 1            | 1        |
| 3           | C     | 1043660-NLA | NLA - Footplate, Aluminum, Deluxe, Black - Left (8" W x 5-3/4" L)                   | 1            | 1        |
| 3           | B     | 44200X025   | Footplate, Aluminum, Extra Large, Black - Right (8" W x 8-1/2" L)                   | 1            | 1        |
| 3           | B     | 44200X026   | Footplate, Aluminum, Extra Large, Black - Left (8" W x 8-1/2" L)                    | 1            | 1        |
| 4           | A     | 1016346     | Cap, End, Tube, Black (7/8")  | 1            | 1        |
| 4           | B     | 1027095     | Plug Button, Black (7/8")   | 1            | 1        |
| 5           | A     | 1006981     | Spring, Tension   | 1            | 1        |
| 5           | C,F   | 1036120-NLA | NLA - Spring  | 1            | 1        |
| 5           | B,I   | 1036123     | Spring, Leaf  | 1            | 1        |
| 6           | G     | 40120X001   | Screw, Phillips Pan Head (#8-32 x 3/8")   | 1            | 2        |
| 7           |       | 1070168     | Plug, Push, Footplate   | 1            | 1        |

NOTE: A - Standard Aluminum (1650)  
 B - Extra Large Aluminum (1350)  
 C - Deluxe Aluminum (1450)  
 D - PH904A  
 E - PHAL4A  
 F - Two per side  
 G - Deluxe takes one  
 H - Part #1023804 is no longer available as a service part individually. See category titled "Front Rigging Reference Chart" and order part #1023924 or 1023925 Pivot/Slide Tube and Footplate Assembly (Deluxe Aluminum - 1450) as a replacement.  
 I - Use one or more leaf springs per assembly as required such that footplate may be moved by hand to the up-right position and remain in place.