


Pivot/Slide Tube and Aluminum Footplate PHW93

Heel loops are not included. See footrest and legrest accessories.
For complete pivot/slide tube and footplate assemblies, see front riggings reference chart.



Pivot/Slide Tube and Aluminum Footplate PHW93

**Heel loops are not included. See footrest and legrest accessories.
For complete pivot/slide tube and footplate assemblies, see front riggings reference chart.**

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1027097	Plug Button, Black (3/4")	1	1
2		1048834	Pivot/Slide Tube w/ Plug Button, Black	1	1
2		1098906	Pivot/Slide Tube w/ Plug Button, Black (2" Longer)	1	1
2		1086482	Pivot/Slide Tube w/ Plug Button, Black (4" Longer)	1	1
3		44200X025	Footplate, Aluminum, Extra Large, Black - Right (8" W x 8-1/2" L)	1	1
3		44200X026	Footplate, Aluminum, Extra Large, Black - Left (8" W x 8-1/2" L)	1	1
4		1027095	Plug Button, Black (7/8")	1	1
5	A	1036123	Spring, Leaf	1	1
6		1025729	Screw, Phillips Pan Head Tap (#8-32 x 3/8")	1	2

NOTE: A - Use one or more leaf springs per assembly as required such that footplate may be moved by hand to the up-right position and remain in place.