





Pivot/Slide Tube w/ Footplate Assy. - Aluminum - A1562/A1564



Small /Medium



Extra Large



Extra-Small

Pivot/Slide Tube w/ Footplate Assy. - Aluminum - A1562/A1564

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		44200X032-NLA	NLA - Footplate, Aluminum, Extra Small, Black (5-3/4" W x 6" L)	1	1
1		44200X059	Footplate, Aluminum, Small, Black (6-3/4" W x 6" L)	1	1
1		44200X027	Footplate, Aluminum, Medium, Black (7-3/4" W x 6" L)	1	1
1		44200X025	Footplate, Aluminum, Extra Large, Black - Right (8" W x 8-1/2" L)	1	1
1		44200X026	Footplate, Aluminum, Extra Large, Black - Left (8" W x 8-1/2" L)	1	1
2	A	1006981	Spring, Tension	1	1
2	C,B	1036123	Spring, Leaf	1	1
3		1025729	Screw, Phillips Pan Head Tap (#8-32 x 3/8")	1	2
4		1048834	Pivot/Slide Tube w/ Plug Button, Black	1	1
5		1027097	Plug Button, Black (3/4")	1	1
6		See Chart	Pivot Slide Tube with Footplate	1	1

NOTE: A - Small and Medium Footplates
 B - Extra Large Footplates
 C - Use one or more leaf springs per assembly as required such that footplate may be moved by hand to the up-right position and remain in place.

#6 Pivot Slide Tube with Footplate

Footplate Size	Chair Width	Standard Length	2" Longer	4" Longer
Extra Small Right	12-13"	1026398-NLA	N/A	N/A
Extra Small Left	12-13"	1026399-NLA	N/A	N/A
Small Right	14-15"	1041542-NLA	N/A	N/A
Small Left	14-15"	1041543-NLA	N/A	N/A
Medium Right	16-20"	1026401	N/A	N/A
Medium Left	16-20"	1026402	N/A	N/A
Extra Large Right	16-20"	1026395	1086515	1086516
Extra Large Left	16-20"	1026396	1086517	1086518