


# Seat Positioning Strap



SCP\_028A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1034625	Strap, Positioning, Seat/Chest w/ Auto Style Push Button Buckle (Eyelet to Eyelet Length 61")	1	1
NOTE: A - Must use existing hardware consisting of 2 ea. Screw, Hex Head Cap (M8 x 1.25mm) and 2 ea. Locknut (M8 x 1.25mm). Hardware is not available as a service part.					