


Miscellaneous Accessories

THE PRODUCTION OF THE P725 IS CURRENTLY SUSPENDED. CONTACT CUSTOMER SERVICE FOR ORDERING OPTIONS.



Miscellaneous Accessories

THE PRODUCTION OF THE P725 IS CURRENTLY SUSPENDED. CONTACT CUSTOMER SERVICE FOR ORDERING OPTIONS.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1035614	Plug Button, Black (5/8")	1	2
2		1027684	Knob, Clamping, Black	1	4
3		1080866	Decal, Oxygen Warning - English	1	1
		1090093-NLA	NLA - Decal, Oxygen Warning - French Canadian	1	1
		1092793-NLA	NLA - Decal, Oxygen Warning - Domestic Spanish	1	1
4		1072004	Pad, Non-Skid Adhesive Back	1	2
5	A	1086780-NLA	NLA - Package, Protective Cap	10	1
6	B	1091306-NLA	NLA - Oxygen Holder Assembly - Deluxe Seats	1	1
7		1088765-NLA	NLA - Fastener Hardware, Accessory Tube Includes: 1 ea. Tighten Knob 1 ea. Screw, Hex Cap (M10 x 25) 1 ea. Nut (M10 x 1.25 x 17)	1	1
	D	2000454-NLA	NLA - Bag, Poly (39" L x 33" W)	1	1
	D	1091305-NLA	NLA - Instruction Sheet, Oxygen Holder	1	1
	D	1083052-NLA	NLA - Carton, Shipper, Oxygen Holder	1	1
	D	P785-NLA	NLA - Basket, Rear - Lynx SX3P/ SX3 and LX3 Plus with Deluxe Seat	1	1
	C,D	P725	Kit, Crutch/Cane Holder	1	1
			Holder, Crutch/Cane - P725	1	1
			Cup, Crutch/Cane Holder, Black	1	1
			Strap, Hook and Loop - P725	1	1
			Plug Button, Black (7/8")	1	1
			Screw w/ Washer (#10 x 7/8")	1	2
	D	P710	Flag, Safety	1	1
	D	1035917-NLA	NLA - Heel Loop (Pair)	2	1
	D	1017621	Heel Loop - Right/Left	1	1
	D,E	1035907	Kit, Heel Loop Attaching Hardware (Pair)	2	1
NOTE: A - Four protective caps are included with the oxygen holder assembly, replacement caps are sold in packages of 10. B - Includes items 1-3, 5 and 6. Fastener hardware part number 1088765 is not included with the oxygen holder assembly. C - THE PRODUCTION OF THE P725 IS CURRENTLY SUSPENDED. CONTACT CUSTOMER SERVICE FOR ORDERING OPTIONS. D - Not Shown E - Consists of: 2 ea. Screw, Phillips Head (1/4-20 x 2-1/4") 2 ea. Spacer (1/4 x 11/32 x 1-3/4") 2 ea. Locknut (1/4-20)					