


# Wheel Locks



SCP\_026A

## Wheel Locks

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1191494-NLA	NLA - Wheel Lock Assembly, Pull to Lock, Black - Right	1	1
1		1191495-NLA	NLA - Wheel Lock Assembly, Pull to Lock, Black - Left	1	1
2		47376X024	Tip, Rubber Wheel Lock Handle, Black	1	1
3		1030949	Screw, Hex Head (1/4-20 x 1-1/2")	1	1
4		1025736	Locknut (1/4-20)	1	1