

## Rear Wheels

**NOTE: If your chair was built with a serial # starting with "03H" or "03I", please contact an Invacare Parts Specialist for replacement parts at 1-800-333-6900 Ext. 3722.**

| Model Number | Wheel Style  | Single Wheel Assembly                      | Tire                                    | Tube/Insert | Bearing   | Bearing Spacer |
|--------------|--|--|---|-------------|-----------|----------------|
| U2222C       | Composite Corded Urethane (24 x 1") 5/8" Axle (9-Spoke)                  | 1121043<br><small>See Note(s) A</small>    | 1118209<br><small>See Note(s) C</small> | N/A         | 40117X006 | 1019695-NLA    |
| 5244BK       | Composite Pneumatic w/Flat Free Insert (24 x 1-3/8") 5/8" Axle (7-Spoke) | 1045281<br><small>See Note(s) B, A</small> | 40176X060                               | 1045280     | 40117X006 | 1019695-NLA    |
| U2244BK      | Composite Urethane (24 x 1-1/4") 5/8" Axle (7-Spoke)                     | 1034130<br><small>See Note(s) A</small>    | 1026904                                 | N/A         | 40117X006 | 1019695-NLA    |

NOTE: A - Includes Tire, Tube/Insert, Bearing and Spacer.  
 B - Required for 450lb weight capacity.  
 C - Corded Tires maybe more difficult to stretch over the wheel rim. The complete wheel assembly that has the tire and wheel rim preassembled maybe advisable

NOTE: Effective 3/27/06, the U2222 wheel assembly is no longer available for this model.

### Rear Wheel and Handrim Assemblies - 24" Wheel with 5/8" Axle

| Wheel Style | Description  | Part Number                             |
|-------------|--|---|
| U2222C      | Assembly, Wheel 24", 5/8" Axle, Composite Urethane - Composite Handrim | 1133298<br><small>See Note(s) A</small> |

NOTE: A - Non-Corded Tire 24 x 1" Low Profile Urethane Tire was replaced with a Corded Tire as of 11/16/09. For chairs manufactured before 11/16/09, it is recommended to order (2) ea. in quantity to match

### Wheel Conversion Kits

| Part Number                             | Description                   | Quantity |
|---|-------------------------------|----------|
| 1121459<br><small>See Note(s) A</small> | Kit, 350Lb Wheel Package (T4) | 1 Pair   |
| 1121460<br><small>See Note(s) B</small> | Kit, 450Lb Wheel Package (T4) | 1 Pair   |

## Rear Wheels

**NOTE: If your chair was built with a serial # starting with "03H" or "03I", please contact an Invacare Parts Specialist for replacement parts at 1-800-333-6900 Ext. 3722.**

NOTE: A - Kit contains 2 - U2222C Wheel assemblies with Alum. Clear Anodized Handrims, 2 - AT903 Caster assemblies and a 350lb weight capacity decal.  
B - Kit contains 2 - 5244BK Wheel assemblies with Alum. Clear Anodized Handrims, 2 - 1228 Caster assemblies and a 450lb weight capacity decal.