

## Non-Projection Handrims and Attaching Hardware

**NOTE: If your chair was built with a serial # starting with "03H" or "03I", please contact an Invacare Parts Specialist for replacement parts at 1-800-333-6900 Ext. 3722.**

As of 8/11/15 Chrome Handrims were replaced by Alum. Clear Anodized and Steel Black 34P was replaced by Alum. Black Anodized Handrims.

### Composite (COM)

Wheel Style	Wheel Type	Handrim Only	Attaching Points	Attaching Hardware Kit
5244BK	24" Heavy Duty Composite Pneumatic w/Flat Free Insert	1064036	7	1065790 See Note(s) A
U2222	24" Composite Urethane Tire	1063919	9	1042074 See Note(s) A
U2222C	24" Composite Urethane Tire	1063919	9	1042074
NOTE: A - Hardware kit contains hardware for (1) handrim. See category titled, "Non-Projection Handrims Attaching Hardware" for more details.				

### Chrome Plated (CHDRM)

Wheel Style	Wheel Type	Handrim Only	Attaching Points	Attaching Hardware Kit
U2222	24" Composite Urethane Tire	1055272-NLA	6	1042074 See Note(s) A
U2222C	24" Composite Urethane Tire	1055272-NLA	6	1042074 See Note(s) A
U2222BK	24" Composite Urethane Tire	1055272-NLA	6	1042074 See Note(s) A
U2244BK	24" Composite Urethane w/Treaded Tire	1054966	7	1025153 See Note(s) A
5244BK	24" Heavy-Duty Composite Pneumatic w/ Flat Free Insert	1054966	7	1025153 See Note(s) A
NOTE: A - Hardware kit contains hardware for (1) handrim. See category titled, "Non-Projection Handrims Attaching Hardware" for more details.				
As of 8/11/15 Chrome Handrims were replaced by Alum. Clear Anodized and Alum. Black 34P was replaced by Alum. Black Anodized Handrims.				

## Non-Projection Handrims and Attaching Hardware

**NOTE: If your chair was built with a serial # starting with "03H" or "03I", please contact an Invacare Parts Specialist for replacement parts at 1-800-333-6900 Ext. 3722.**

As of 8/11/15 Chrome Handrims were replaced by Alum. Clear Anodized and Steel Black 34P was replaced by Alum. Black Anodized Handrims.

### Aluminum - Clear Anodized

Wheel Style	Wheel Type	Handrim Only	Attaching Points	Attaching Hardware Kit
U2222	24" Composite Urethane Tire	1183681	6	1042074 See Note(s) A
U2222C	24" Composite Urethane Tire	1183681	6	1042074 See Note(s) A
U2222BK	24" Composite Urethane Tire	1183681	6	1042074 See Note(s) A
U2244BK	24" Composite Urethane w/Treaded Tire	1183680	7	1025153 See Note(s) A
5244BK	24" Heavy-Duty Composite Pneumatic w/ Flat Free Insert	1183680	7	1025153 See Note(s) A

NOTE: A - Hardware kit contains hardware for (1) handrim. See category titled, "Non-Projection Handrims Attaching Hardware" for more details.

### Steel - Flat Black 34P

Wheel Type	Handrim Only	Attaching Points	Attaching Hardware Kit
U2222 24" Composite Urethane Tire	1120063-NLA See Note(s) B	6	1042074 See Note(s) A
U2222C	1120063-NLA See Note(s) B	6	1042074 See Note(s) A
U2244BK 24" Composite Urethane Treaded Tire	1119821-NLA See Note(s) C	7	1025153 See Note(s) A
5244BK 24" Heavy Duty Composite w/Flat Free Insert	1119821-NLA See Note(s) C	7	1025153 See Note(s) A

NOTE: A - Hardware kit contains enough hardware for (1) handrim. See category titled, "Non-Projection Handrims Attaching Hardware" for more details.  
 B - NLA. See Black Anodized Handrim Chart and order part #1186220.  
 C - NLA. See Black Anodized Handrim Chart and order part #1186219.

NOTE: When replacing handrims for chairs built before 12/17/99 order in pairs. Handrim color has changed from Pewter to Flat Black.  
 NOTE: As of 8/11/15 Chrome Handrims were replaced by Alum. Clear Anodized and Alum. Black 34P was replaced by Alum. Black Anodized Handrims.

## Non-Projection Handrims and Attaching Hardware

**NOTE: If your chair was built with a serial # starting with "03H" or "03I", please contact an Invacare Parts Specialist for replacement parts at 1-800-333-6900 Ext. 3722.**

As of 8/11/15 Chrome Handrims were replaced by Alum. Clear Anodized and Steel Black 34P was replaced by Alum. Black Anodized Handrims.

### Aluminum - Black Anodized

Wheel Type	Handrim Only	Attaching Points	Attaching Hardware Kit
U2222 24" Composite Urethane Tire	1186220	6	1042074 See Note(s) A
U2222C	1186220	6	1042074 See Note(s) A
U2244BK 24" Composite Urethane Treaded Tire	1186219	7	1025153 See Note(s) A
5244BK 24" Heavy Duty Composite w/Flat Free Insert	1186219	7	1025153 See Note(s) A
NOTE: A - Hardware kit contains enough hardware for (1) handrim. See category titled, "Non-Projection Handrims Attaching Hardware" for more details.			