




Seat Upholstery



R32 AND R51 (Before 11/15/02)
Seat Upholstery - Vinyl



R32 and R51 (After 11/14/02)
Seat Upholstery - Nylon



R32LX/R50LX/R51LX/R51LXP
Seat Upholstery - Nylon

Seat Upholstery

Lightweight Vinyl Seat Upholstery

R32 and R51 (Before 11/15/02)

Chair Width	16" Depth	17" Depth	18" Depth
14"	8881018200-NLA	8881052296	8881052301-NLA
16"	8881016920-NLA	8881020690-NLA	8881023360-NLA
18"	8881016918-NLA	8881020691-NLA	8881023359-NLA
20"	8881019577-NLA	8881052300-NLA	8881052307-NLA

Note: When ordering, please specify U05 Black for Lightweight/Vinyl or U240 Black for Nylon Upholstery.

Nylon Seat Upholstery

R32 and R51 (After 11/14/02)

Chair Width	16" Depth	17" Depth	18" Depth
14"	8881057199-NLA	8881057206-NLA	8881057213-NLA
16"	1016321 See Note(s) A, A1	1016322-NLA See Note(s) A, A1	1016323-NLA See Note(s) A1, A
18"	1016327 See Note(s) A, A1	1016328 See Note(s) A, A1	1016329 See Note(s) A, A1
20"	1017643 See Note(s) A, A1	1017646-NLA See Note(s) A1, A	1017649 See Note(s) A, A1

NOTE: A - After 2/2/16 Crossbrace and Seat Upholstery Holes were revised to a larger size. One hole was added per side and attaching hardware Screw changed from (#10 x 7/8" Screw w/ Washer) to (#12 x 7/8" Screw).
A1 - Wheelchairs manufactured Before 2/3/16 will require new Crossbraces, Seat Upholstery and (Hardware After 2/2/16) to be ordered.

Note: When ordering, please specify U05 Black for Lightweight/Vinyl and U240 Black for Nylon Upholstery.

Seat Upholstery

Nylon Seat Upholstery

R32LX/R50LX/R51LX/R51LXP

Chair Width	16" Depth	17" Depth	18" Depth
14"	1016315 See Note(s) A, A1	1016316-NLA See Note(s) A1, A	1016317-NLA See Note(s) A, A1
16"	1016321 See Note(s) A, A1	1016322-NLA See Note(s) A, A1	1016323-NLA See Note(s) A, A1
18"	1016327 See Note(s) A1, A	1016328 See Note(s) A1, A	1016329 See Note(s) A1, A
20"	1017643 See Note(s) A, A1	1017646-NLA See Note(s) A1, A	1017649 See Note(s) A1, A

NOTE: A - After 2/2/16 Crossbrace and Seat Upholstery Holes were revised to a larger size. One hole was added per side and attaching hardware Screw changed from (#10 x 7/8" Screw w/ Washer) to (#12 x 7/8" Screw).
A1 - Wheelchairs manufactured Before 2/3/16 will require new Crossbraces, Seat Upholstery and (Hardware After 2/2/16) to be ordered.

Note: When a different material or color is needed, please see special SM-060831091250.

Hardware for Seat Upholstery

Part Number	Description	16" Depth Quantity Required	17"-18" Depth Quantity Required
1076284 See Note(s) A	Single Pack Screw with Washer (#10 x 7/8") Used Before 2/3/16	8	10
1057950 See Note(s) A	Multi Pack Screws with Washers (#10 x 7/8") Used Before 2/3/16	1 Pkg of 8	N/A
1050393-NLA See Note(s) A	Multi Pack Screws with Washers (#10 x 7/8") Used Before 2/3/16	N/A	1 Pkg of 10
1054978 See Note(s) C, C1, B	Multi Pack Screws and Washers (#12 x 7/8") Used After 2/2/16	1 Pkg of 10	N/A
1131378 See Note(s) C1, B, C	Multi Pack Screws and Washers (#12 x 7/8") Used After 2/2/16	N/A	1 Pkg of 12

NOTE: A - Used Before 2/3/16.
B - Used After 2/2/16.
C - After 2/2/16 Crossbrace and Seat Upholstery Holes were revised to a larger size. One hole was added per side and attaching hardware Screw changed from (#10 x 7/8" Screw w/ Washer) to (#12 x 7/8" Screw).
C1 - Wheelchairs manufactured Before 2/3/16 will require new Crossbraces, Seat Upholstery and (Hardware After 2/2/16) to be ordered.