


Mattress Support, Leg End



SBP_068A

Mattress Support, Leg End

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,B	SP1584447	Mattress Support, Complete, Leg End	1	1
1	A,C	SP1587786	Mattress Support, Complete, Leg End	1	1
2		SP1435195	Kit, Rivet Replacement	1	1
3		SP1417418	Kit, Ratchet and Pin, Four Sections	1	1
4	B	SP1423882	Top Frame, Leg Section	1	1
4	C	SP1587787	Top Frame, Leg Section	1	1
5	B	SP319114	Leg Section	1	1
5	C	SP1587788	Leg Section	1	1
6		SP217996-7035	Handle for Leg Section	1	1
7		SP1574230	Mattress Holder (Pair)	2	1
8		SP1584437	Screw, Finger (M8 x 13) (Pair)	2	1

NOTE: A - Includes items 2 - 8.
 B - Before 4/1/16
 C - After 3/31/16