


# Front Riggings - 70 Degree Hanger (70SA)



SCM\_129B

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	#1. 70 Degree Hanger - Complete	1	2
2		SP1584383	Kit, Hanger Locking Nose	1	2
3		SP1584381	Kit, Hanger Pivot	1	2
4		SP1584380	Kit, Hanger Locked	1	2
5		SP1587483	Kit, Hanger Stopper	1	2
6		SP1584382	Calf Strap Holder and Ball Screw	1	2
7		SP1424571	Screw, Stop (M10 x 23)	1	2

NOTE: A - Includes items 2, 3, 4, 5, 6, and 7.

## #1. 70 Degree Hanger - Complete

Color	Part Number	Color	Part Number
Pearl White	SP1584371-0637	Hypnotic Purple	SP1584371-0645
Anthracite Black	SP1584371-0609	Pearl Grey	SP1584371-0515
Pop Blue	SP1584371-0643	Matte Black	SP1584371-0646
Deep Red	SP1584371-0607	Intense Yellow	SP1584371-0647
Extreme Green	SP1584371-0644	Vibrant Orange	SP1584371-0648