


Wheels - Rear - Supports



SCM_161B

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	SP1554604-0005	Kit, Wheel Assembly, Rear Support w/o Accent - Left	1	1
1	A	SP1555921-0005	Kit, Wheel Assembly, Rear Support w/o Accent - Right	1	1
1	A	SP1554604-0641	Kit, Wheel Assembly, Rear Support, Blue Accent - Left	1	1
1	A	SP1555921-0641	Kit, Wheel Assembly, Rear Support, Blue Accent - Right	1	1
1	A	SP1554604-0642	Kit, Wheel Assembly, Rear Support, Purple Accent - Left	1	1
1	A	SP1555921-0642	Kit, Wheel Assembly, Rear Support, Purple Accent - Right	1	1
2		SP1601541	Kit, Hardware Wheel Assembly, Rear Support	1	1
3		SP1601542	Kit, Wheel Assembly, Rear Bearing	1	1
4		SP1584797	Kit, Camber, 0 - 2 Degree (Pair)	1	1
5		SP1584798	Kit, Camber, 3 - 4 Degree (Pair)	1	1

NOTE: A - Includes items 2, 3, 4, and 5.