


Wheel Locks - Attendant Operated, Hub Brake



SCM_146A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		SP1601543	Kit, Attendant Hub Brake, Complete (w/o Wheels)	1	1
2		SP1601544	Kit, Hub Brake Plate w/ Spacer - Left	1	1
2		SP1601545	Kit, Hub Brake Plate w/ Spacer - Right	1	1
3		SPM30071	Cable Clamp Bolt w/ Nut and Washer	1	1
4		SP1425473	Handle, Back Cane (7/8)	1	1
5		SP1425530	Insert, Push Handle	1	1
6		SP1428603	Screw Tfb (4,8 x 19 A2r Tx25 DIN7983)	1	1
7		SP5221363	Hub Brake Adaptation, Back Cane, Fixed	1	1
8		SP1420527	Lever, Brake, Black (7/8 Dia.)	1	1
9		SP1564821	Attendant Hub Brakes, Complete Mechanism (Pair)	1	1