


# Back Cane Attaching Hardware

## Fixed Back (FXBACK)



SCM\_197A

## Back Cane Attaching Hardware Fixed Back (FXBACK)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1164103	Kit, Back Cane Attaching Hardware, Fixed	1	1
2			Plate, Fixed Back, Black	1	2
3			Screw, Hex Head (1/4-28 x 2")	1	1
4			Bracket, Backrest, Black	1	2
5			Screw, Hex Head (1/4-28 x 2")	1	1
6			Screw, Hex Head (1/4-28 x 2-1/4")	1	2
7			NLA - Washer, Steel (1/4 x 1/2 x 3/64")	1	2
8			Uni-Washer, Dual	1	1
9			Washer, Steel (1/4 x 5/8 x 3/32")	1	2
10			Nut, Hex Nylon (1/4-28)	1	1
11			Locknut (1/4-28)	1	1
12			Locknut, Thin (1/4-28)	1	1
13			Uni-Washer, Coved	1	2
14			NLA - Bushing, Nylon (1/4 x 3/8 x 1-1/4")	1	1
15			Screw, Socket Button Head (1/4-20 x 1-3/4")	1	1
16			Locknut, Thin (1/4-20)	1	1
17			Nut, Hex Nylon (1/4-28)	1	1

NOTE: A - Includes items 2 - 17. Hardware quantity accommodates one side only.