


# Back Cane Attaching Hardware

## Fold Down Back (FLDBACK)



## Back Cane Attaching Hardware Fold Down Back (FLDBACK)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1164102	Kit, Back Cane Attaching Hardware, Fold Down	1	1
2			Bracket, Latch Back	1	1
3			NLA - Spring, Compression	1	1
4			Spacer, Adjuster	1	1
5			Plate, Adjuster, Backrest	1	2
6			Uni-Washer, Coved	1	1
7			Screw, Socket Button Head (1/4-20 x 1-3/4")	1	1
8			Bracket, Backrest, Black	1	2
9			Screw, Hex Head (1/4-28 x 2")	1	1
10			Screw, Hex Head (1/4-28 x 2")	1	1
11			Screw, Hex Head (1/4-28 x 2-1/4")	1	2
12			NLA - Washer, Steel (1/4 x 1/2 x 3/64")	1	2
13			Uni-Washer, Dual	1	1
14			Washer (1/4 x 5/8 x 3/32")	1	2
15			Nut, Hex Nylon (1/4-28)	1	1
16			Locknut (1/4-28)	1	1
17			Locknut, Thin (1/4-28)	1	1
18			Locknut, Thin (1/4-20)	1	2
19			NLA - Spacer (1/4 x 3/8 x 3/16")	1	2
20			Screw, Socket Button Head (1/4-20 x 1-5/8")	1	1
21			Nut, Hex Nylon (1/4-28)	1	1
NOTE: A - Includes items 2 - 21. Hardware quantity accommodates one side only.					