


Arms - Non-Flip Back Attaching Hardware (HARMHT Option)


HARMHT - High Mount Attaching Hardware Setting Left/Right Rear Arm Socket



<p>Right Side View Shown</p> <p>5</p> <p>Position A</p>	SDRI216	D15
	SDRI519	D15
	SDRI822	D18

<p>Right Side View Shown</p> <p>5</p> <p>Position B</p>	SDRI216	D16
	SDRI519	N/A
	SDRI822	N/A

HARMHT - High Mount Attaching Hardware Setting Right/Left Rear Arm Socket



<p>Right Side View Shown</p> <p>3</p> <p>Position A</p>	SDRI216	N/A
	SDRI519	D16
	SDRI822	D19

<p>Right Side View Shown</p> <p>3</p> <p>Position B</p>	SDRI216	N/A
	SDRI519	D17 D18 D19
	SDRI822	D20 D21 D22

Arms - Non-Flip Back Attaching Hardware (HARMHT Option)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,C	1164104	Kit, Arm Attaching Hardware - Right	1	1
2			Socket, Arm, Rear, Black - Right	1	1
3			Socket, Arm, Front, Black - Right	1	1
6			NLA - Screw (1/4-20 x 1-3/4")	1	1
7			Nut, Hex with Nylon Cap (1/4-20)	1	1
8	B,C	1164105	Kit, Arm Attaching Hardware - Left	1	1
4			Socket, Arm, Rear, Black - Left	1	1
5			Socket, Arm, Front, Black - Left	1	1
6			NLA - Screw (1/4-20 x 1-3/4")	1	1
7			Nut, Hex with Nylon Cap (1/4-20)	1	1
		See Chart	Arm Settings and Hardware (HARMHT Option)	1	1

NOTE: A - Kit includes items 2, 3, 6, and item 7. Parts are not sold separately.
 B - Kit includes items 4, 5, 6, and item 7. Parts are not sold separately.
 C - For further information, see "Arm Settings and Hardware (HARMHT Option)" Chart.

Arm Settings and Hardware (HARMHT Option)

See Illustration for Corresponding Information

HARMHT Option	Rear Arm Socket Left (Item 5)	Seat Depth (In.)	Arm Option
1164105	SDR1216	D15	
Left	SDR1519	D15	
Position A	SDR1822	D18	
1164105	SDR1216	D16	
Left	SDR1519	N/A	
Position B	SDR1822	N/A	
	Item 5 Mounted on Right Side		FXDESK - FXFULL ADJDESK - ADJFULL
HARMHT Option	Rear Arm Socket Right (Item 3)	Seat Depth (In.)	Arm Option
1164104	SDR1216	N/A	
Right	SDR1519	D16	
Position A	SDR1822	D19	
1164104	SDR1216	N/A	
Right	SDR1519	D17, D18, D19	
Position B	SDR1822	D20, D21, D22	
	Item 3 Mounted on Right Side		FXDESK - FXFULL ADJDESK - ADJFULL