


Frame, Wire (Grid) Deck



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1124075-NLA	NLA - Frame, Bed & Deck, Taupe - Wire	1	1
2		026507	Pan Deck - Seat Section Assembly	1	1
3		026623	Rivet (1/4")	1	6
4		1125002	Cap, Angle - Left (2 x 1-1/2 x 1/8")	1	1
5		1125015	Cap, Angle - Right (2 x 1-1/2 x 1/8")	1	1
6	A	1124072-NLA	NLA - Frame & Wire Deck Assembly, Taupe	1	1
	B	1124029-NLA	NLA - Kit, Corner Bracket, Taupe	1	1
	B	1124031	Kit, Mattress Keeper Includes: 4 ea. Mattress Keepers, 8 ea. Nuts, 8 ea. Screws, and 1 ea. Poly Bag	1	1

NOTE: A - Includes: Items 1-5 and Items Not Shown.
B - Not Shown.