



Carriage Floor Lock



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1131077	Foot, Rubber	1	4
2		1012756	Plug Button, Black (1-1/4")	1	4
3		020717	Glide, Cap	1	2
4	A	1124074	Carriage Floor Lock Assembly	1	1
5		009077	Kit, Caster	1	1
6		026111	Fixed Link Assembly, Knee	1	1

NOTE: A - Includes: Items 1-3 and Frame shown.