


Manual Pull Tube and Crank Handle



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1124064	Manual, Head Crank and Pull Tube Assembly	1	1
2		1124063	Manual, Hi-Lo Crank and Pull Tube Assembly	1	1
3		1124061	Manual, Foot Crank and Pull Tube Assembly	1	1
4		001008	Washer, Split Lock (5/16")	1	6
5		001741	Screw, Hex Head (5/16-18 x 3/4")	1	6
6		005952	Pin, Clevis (5/16 x 2-1/8")	1	3
7		51000M743	Pin, Hitch (3/8 x 1-3/16")	1	3
8		51000M191	Cover, Pull Tube (Head/Foot)	1	3
9		001405	Bushing, Flanged, Nylon	1	3
10		001443	Bushing, Shoulder - Pivot	1	3
11	A	1124077	Kit, Manual 3 Crank	1	1

NOTE: A - Includes: Items 1-10.