


Frame, Wire (Grid) Deck 820DLX, 820-3MDLX & SC900DLX



Frame, Wire (Grid) Deck 820DLX, 820-3MDLX & SC900DLX

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1137371	Frame and Wire Deck w/ Seat Pan	1	1
2		1125015	Cap, Angle - Right (2 x 1-1/2 x 1/8")	1	1
3		1125002	Cap, Angle - Left (2 x 1-1/2 x 1/8")	1	1

NOTE: A - Includes: Items 2-3