



Seat Upholstery and Hardware Kits

NOTE: If your chair was built with a serial # starting with "03H" or "03I", please contact an Invacare Parts Specialist for replacement parts at 1-800-333-6900 Ext. 3722.

NOTE: After 10/11/16 additional hole was added to the rear corners of the seat upholstery and accommodates all wheelchairs manufactured before 10/12/16 as well. All Crossbrace Seat Rails in the field already have the additional upholstery hole and do not require replacement of Crossbraces.



Seat Upholstery - Nylon



Seat Upholstery - Vinyl

Seat Upholstery and Hardware Kits

NOTE: If your chair was built with a serial # starting with "03H" or "03I", please contact an Invacare Parts Specialist for replacement parts at 1-800-333-6900 Ext. 3722.

NOTE: After 10/11/16 additional hole was added to the rear corners of the seat upholstery and accommodates all wheelchairs manufactured before 10/12/16 as well. All Crossbrace Seat Rails in the field already have the additional upholstery hole and do not require replacement of Crossbraces.

Chair Width x Depth	Vinyl Upholstery	Nylon Upholstery
18" x 18"	8881127512 See Note(s) A, B	1127579 See Note(s) B
18" x 20"	8881127588 See Note(s) A, B	1127580 See Note(s) B
20" x 18"	8881135515 See Note(s) B, A	1127536 See Note(s) B
20" x 20"	8881127590 See Note(s) A, B	1127582 See Note(s) B
22" x 18"	8881127591 See Note(s) B, A	1127538 See Note(s) B
22" x 20"	8881127592 See Note(s) B, A	1127584 See Note(s) B
24" x 18"	8881127593 See Note(s) B, A	1127585 See Note(s) B
24" x 20"	8881127594 See Note(s) B, A	1127586 See Note(s) B

NOTE: A - Specify upholstery color when ordering. See category titled, "Upholstery and Frame Colors."
B - Upholstery kits include hardware. 18" deep chairs have six mounting holes and 20" deep chairs have seven mounting holes.

Hardware - Seat Upholstery

Part Number	Description	Quantity
1194392 See Note(s) A	14 ea. Screw, Phillips Oval Head (#12-14 x 7/8") 14 ea. Washer (1/4 x 3/4 x 3/32")	1 package of 14 Screws and 14 Washers

NOTE: A - Hardware is not sold individually.