

# Frames without Hardware

## 9000SL

Chair Style	Side	Chrome Part Number	Flat Black - 34P Part Number	Color Part Number
Adult/Hemi Before 7/12/99	Right	1056598-NLA	1066267-NLA	8881066269-NLA See Note(s) A, C
Adult/Hemi Before 7/12/99	Left	1056599-NLA	1066268-NLA	8881066270-NLA See Note(s) A, C
Adult/Hemi After 7/11/99	Right	1091096-NLA See Note(s) C1, C	1091098 See Note(s) C1, C	8881091100 See Note(s) A, C1, C
Adult/Hemi After 7/11/99	Left	1091097-NLA See Note(s) C1, C	1091099 See Note(s) C1, C	8881091101 See Note(s) A, C1, C
Adjustable Angle Back	Right	1063910-NLA See Note(s) C1, C, B	1064774 See Note(s) C1, C, B	8881063917 See Note(s) A, C1, B, C
Adjustable Angle Back	Left	1063911-NLA See Note(s) C1, C, B	1064775 See Note(s) C1, B, C	8881063918 See Note(s) A, C1, C, B

NOTE: A - Specify frame finish when ordering, see category titled, "Frame Finishes, Upholstery Colors, Decals and Packaging."  
 B - Adjustable angle back frame is only available with adjustable height conventional style arms. Adjustable angle back frames are the same regardless of the chair's seat-to-floor height. All frames include caster headtube and bearings.  
 C - Hinge Pin Replacement Hardware Kit Located on Frame Hardware page is used on Side Frames w/ welded-on Hinge Pins manufactured before 4/12/13. Side Frames after 4/11/13 have Hinge Pin Saddle Plates w/ square holes to attach REMOVABLE Non-Welded Hinge Pins.  
 C1 - All Side Frames ordered after 4/11/13 require new Side Frame Hardware that will include the removable Hinge Pins also. Removable Hinge Pins are not sold separately.

NOTE: When ordering frames for chairs built before 10/3/97, both left and right side frames must be ordered due to a style change.

## 9000XT

Chair Style	Side	Chrome Part Number	Flat Black - 34P Part Number Before 10/12/16	Flat Black - 34P Part Number After 10/11/16	Color Part Number Before 10/12/16	Color Part Number After 10/11/16
Adult/Hemi/ Super Hemi	Right	1056981-NLA See Note(s) B1, B	1056977 See Note(s) B1, B	1194284 See Note(s) B1, B	8881056979 See Note(s) A, B1, B	8881194283 See Note(s) A, B, B1
Adult/Hemi/ Super Hemi	Left	1056982-NLA See Note(s) B1, B	1056978 See Note(s) B1, B	1194276 See Note(s) B, B1	8881056980 See Note(s) B1, B, A	8881194279 See Note(s) B1, B, A
Adult/Hemi/ Super Hemi with Optional Headtube	Right	1091106-NLA See Note(s) B1, B	1091108 See Note(s) B1, B	1091108 See Note(s) B1, B	8881091070 See Note(s) B1, B, A	N/A
Adult/Hemi/ Super Hemi with Optional Headtube	Left	1091107-NLA See Note(s) B1, B	1091109 See Note(s) B1, B	1091109 See Note(s) B1, B	8881091071 See Note(s) B1, B, A	NA
Super Low 16-17" Deep	Right	1063167-NLA See Note(s) B1, B	1063163 See Note(s) B1, B	1063163 See Note(s) B1, B	8881063165 See Note(s) B1, B, A	8881194263 See Note(s) B1, A, B
Super Low 16-17" Deep	Left	1063168-NLA See Note(s) B1, B	1063164 See Note(s) B1, B	1063164 See Note(s) B1, B	8881063166 See Note(s) B, A, B1	8881194260 See Note(s) A, B1, B
Super Low 18" Deep	Right	1079219-NLA See Note(s) B1, B	1079221 See Note(s) B, B1	1079221 See Note(s) B1, B	8881079223 See Note(s) B1, B, A	8881194302 See Note(s) B1, B, A
Super Low 18" Deep	Left	1079220-NLA See Note(s) B1, B	1079222 See Note(s) B1, B	1079222 See Note(s) B1, B	8881079224 See Note(s) B1, B, A	8881194290 See Note(s) B1, B, A
Adjustable Angle Back	Right	1063910-NLA See Note(s) B1, B	1064774 See Note(s) B1, B	1064774 See Note(s) B1, B	8881063917 See Note(s) B1, B, A	8881194259 See Note(s) B1, B, A

# Frames without Hardware

9000XT

Chair Style	Side	Chrome Part Number	Flat Black - 34P Part Number Before 10/12/16	Flat Black - 34P Part Number After 10/11/16	Color Part Number Before 10/12/16	Color Part Number After 10/11/16
Adjustable Angle Back	Left	1063911-NLA See Note(s) B1, B	1064775 See Note(s) B1, B	1064775 See Note(s) B1, B	8881063918 See Note(s) B1, A, B	8881194257 See Note(s) B, B1, A
<p>NOTE: A - Specify frame finish when ordering, see category titled, "Frame Finishes, Upholstery Colors, Decals and Packaging."            B - Hinge Pin Replacement Hardware Kit Located on Frame Hardware page is used on Side Frames w/ welded-on Hinge Pins manufactured before 4/12/13. Side Frames after 4/11/13 have Hinge Pin Saddle Plates w/ square holes to attach REMOVABLE Non-Welded Hinge Pins.            B1 - All Side Frames ordered after 4/11/13 require new Side Frame Hardware that will include the removable Hinge Pins also. Removable Hinge Pins are not sold separately.</p> <p>NOTE: When ordering frames for chairs built before 10/3/97, both left and right side frames must be ordered due to a style change.</p>						

# Frames without Hardware

## 9000XT Recliner

Chair Style	Side	Chrome Part Number	Flat Black - 34P Part Number	Color Part Number Before 10/12/16	Color Part Number After 10/11/16
With Out Lift Kit	Right	1072058-NLA See Note(s) B, B1	1072054 See Note(s) B, B1	8881072056 See Note(s) B, A, B1	8881194341 See Note(s) B1, A, B
With Out Lift Kit	Left	1072059-NLA See Note(s) B, B1	1072055 See Note(s) B, B1	8881072057 See Note(s) B1, B, A	8881194334 See Note(s) B1, A, B
With Lift Kit	Right	1074637-NLA	1074633 See Note(s) B1, B	8881074635 See Note(s) B, A, B1	8881194342 See Note(s) B1, A, B
With Lift Kit	Left	1074638-NLA See Note(s) B, B1	1074634 See Note(s) B, B1	8881074636 See Note(s) B, A, B1	8881194305 See Note(s) B, A, B1

NOTE: A - Specify frame finish when ordering, see category titled, "Frame Finishes, Upholstery Colors, Decals and Packaging."  
 B - Hinge Pin Replacement Hardware Kit Located on Frame Hardware page is only used on Side Frames w/ welded-on Hinge Pins manufactured before 4/12/13. Side Frames after 4/11/13 have Hinge Pin Saddle Plates w/ square holes to attach REMOVABLE Non-Welded Hinge  
 B1 - All Side Frames ordered after 4/11/13 require new Side Frame Hardware that will include the removable Hinge Pins also. Removable Hinge Pins are not sold separately.

NOTE: When ordering frames for chairs built before 10/3/97, both left and right side frames must be ordered due to a style change.

## 9000XDT

Chair Style	Side	Chrome Part Number	Flat Black - 34P Part Number	Color Part Number
Adult/Hemi 16-18" Deep	Right	1074886-NLA See Note(s) B1, B	1074888 See Note(s) B, B1	8881074890 See Note(s) B, A, B1
Adult/Hemi 16-18" Deep	Left	1074887-NLA See Note(s) B, B1	1074889 See Note(s) B1, B	8881074891 See Note(s) A, B, B1
Adult/Hemi 19-20" Deep	Right	1074892-NLA See Note(s) B1, B	1074894 See Note(s) B1, B	8881074896 See Note(s) B1, B, A
Adult/Hemi 19-20" Deep	Left	1074893-NLA See Note(s) B1, B	1074895 See Note(s) B1, B	8881074897 See Note(s) B1, B, A
Tall 16-18" Deep	Right	1074898-NLA See Note(s) B1, B	1074900 See Note(s) B1, B	8881074902 See Note(s) B1, B, A
Tall 16-18" Deep	Left	1074899-NLA See Note(s) B1, B	1074901 See Note(s) B1, B	8881074903 See Note(s) B1, B, A
Tall 19-20" Deep	Right	1074904-NLA See Note(s) B1, B	1074906 See Note(s) B1, B	8881074908 See Note(s) B, A, B1
Tall 19-20" Deep	Left	1074905-NLA See Note(s) B1, B	1074907 See Note(s) B1, B	8881074909 See Note(s) B1, B, A

NOTE: A - Specify frame finish when ordering, see category titled, "Frame Finishes, Upholstery Colors, Decals and Packaging."  
 B - Hinge Pin Replacement Hardware Kit Located on Frame Hardware page is only used on Side Frames w/ welded-on Hinge Pins manufactured before 4/12/13. Side Frames after 4/11/13 have Hinge Pin Saddle Plates w/ square holes to attach REMOVABLE Non-Welded Hinge  
 B1 - All Side Frames ordered after 4/11/13 require new Side Frame Hardware that will include the removable Hinge Pins also. Removable Hinge Pins are not sold separately.

NOTE: When ordering frames for chairs built before 10/3/97, both left and right side frames must be ordered due to a style change.

NOTE: When ordering replacement frames for a previous special frame modification, please contact an Invacare Parts Specialist at 1-800-333-6900 extension 3722.