


Front Wheels - Forks - Aluminum and Titanium Frame



Front Wheels - Forks - Aluminum and Titanium Frame

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1562560-0004	Supporter, Front Wheel (Pair), Black	2	1
2		1562561-0004	Supporter, Front Wheel - Left, Black	1	1
3		1562562-0004	Supporter, Front Wheel - Right, Black	1	1
4		SP1590122	Kit, Fork, Front Pin Bent (Pair)	2	1
5		SP1590123	Kit, Fork, Front Pin, Straight (Pair)	2	1
6		See Chart	Front Wheels - Fork - Aluminum and Titanium Frame	2	1
7		See Chart	Front Wheels - Fork - Aluminum and Titanium Frame	2	1
8		See Chart	Front Wheels - Fork - Aluminum and Titanium Frame	2	1
9		1510219	Kit, Bearing	1	1
		See Chart	Configuration Table, Front Fork Mounted, 75 Degree Frame - Short	1	1
		See Chart	Configuration Table, Front Fork Mounted, 75 Degree Frame - Long	1	1
		See Chart	Configuration Table, Front Fork Mounted, 85 Degree Frame - Short	1	1
		See Chart	Configuration Table, Front Fork Mounted, 85 Degree Frame - Long	1	1

Front Wheels - Fork - Aluminum and Titanium Frame

Color	Item #6 (Small)	Item #7 (Medium)	Item #8 (Large)
Chrome	1562565-0001	1562566-0001	1562567-0001
Black Anodized	1562565-0004	1562566-0004	1562567-0004
Champagne	1562565-0182	1562566-0182	1562567-0182
Techno-Anthrazit	1562565-0400	1562566-0400	1562567-0400
B - Blue	1562565-0404	1562566-0404	1562567-0404
Comet Red	1562565-0440	1562566-0440	1562567-0440
Arctic Blue	1562565-0441	1562566-0441	1562567-0441

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 75 Degree Frame - Short

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
450	390	L	M	M	M	S	S
450	400	M	M	S	M	S	S
450	410	M	M	S	M	S	S
450	420	M	S	S	M	S	S
450	430	M	S	S	M	S	S
450	440	M	S	S	M	S	S
450	450	S	S	S	S	S	S
460	390	L	L	M	L	M	M
460	400	L	L	M	M	M	S
460	410	L	M	M	M	M	S
460	420	L	M	M	M	M	S
460	430	M	M	S	M	M	S
460	440	M	M	S	M	S	S
460	450	M	M	S	M	S	S
460	460	M	S	S	M	S	S
470	400	L	L	L	L	M	M
470	410	L	L	M	L	M	M
470	420	L	L	M	L	M	M
470	430	L	M	M	L	M	M
470	440	L	M	M	M	M	S
470	450	L	M	M	M	M	S
470	460	M	M	S	M	M	S
470	470	M	M	S	M	M	S
480	410	L	L	L	L	L	M
480	420	L	L	L	L	L	M
480	430	L	L	L	L	L	M
480	440	L	L	M	L	M	M
480	450	L	L	M	L	M	M
480	460	L	L	M	L	M	M
480	470	L	M	M	L	M	M
480	480	L	M	M	L	M	M
490	420	-	L	L	L	L	L
490	430	-	L	L	L	L	L
490	440	-	L	L	L	L	M
490	450	-	L	L	L	L	M
490	460	-	L	L	L	L	M
490	470	L	L	M	L	L	M
490	480	L	L	M	L	L	M
490	490	L	L	M	L	L	M
500	430	-	-	L	-	L	L
500	440	-	-	L	-	L	L
500	450	-	-	L	-	L	L
500	460	-	-	L	-	L	L
500	470	-	L	L	-	L	L

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 75 Degree Frame - Short

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
500	480	-	L	L	-	L	L
500	490	-	L	L	-	L	L
510	440	-	-	-	-	-	L
510	450	-	-	-	-	-	L
510	460	-	-	-	-	-	L
510	470	-	-	-	-	-	L
510	480	-	-	L	-	-	L
510	490	-	-	L	-	-	L
520	450	-	-	-	-	-	-
520	460	-	-	-	-	-	-
520	470	-	-	-	-	-	-
520	480	-	-	-	-	-	-
520	490	-	-	-	-	-	-
530	460	-	-	-	-	-	-
530	470	-	-	-	-	-	-
530	480	-	-	-	-	-	-
530	490	-	-	-	-	-	-
540	470	-	-	-	-	-	-
540	480	-	-	-	-	-	-
540	490	-	-	-	-	-	-

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 75 Degree Frame - Long

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
450	390	S	S	S	S	S	S
450	400	S	S	S	S	S	S
450	410	S	S	S	S	S	S
450	420	S	S	S	S	S	S
450	430	S	S	S	S	S	S
450	440	S	S	S	S	S	S
450	450	S	S	S	S	S	S
460	390	M	S	S	S	S	S
460	400	S	S	S	S	S	S
460	410	S	S	S	S	S	S
460	420	S	S	S	S	S	S
460	430	S	S	S	S	S	S
460	440	S	S	S	S	S	S
460	450	S	S	S	S	S	S
460	460	S	S	S	S	S	S
470	400	M	S	S	S	S	S
470	410	M	S	S	S	S	S
470	420	S	S	S	S	S	S
470	430	S	S	S	S	S	S
470	440	S	S	S	S	S	S
470	450	S	S	S	S	S	S
470	460	S	S	S	S	S	S
470	470	S	S	S	S	S	S
480	410	M	M	S	M	S	S
480	420	M	M	S	S	S	S
480	430	M	S	S	S	S	S
480	440	M	S	S	S	S	S
480	450	S	S	S	S	S	S
480	460	S	S	S	S	S	S
480	470	S	S	S	S	S	S
480	480	S	S	S	S	S	S
490	420	L	M	M	M	S	S
490	430	M	M	S	M	S	S
490	440	M	M	S	M	S	S
490	450	M	S	S	M	S	S
490	460	M	S	S	S	S	S
490	470	M	S	S	S	S	S
490	480	S	S	S	S	S	S
490	490	S	S	S	S	S	S
500	430	L	L	M	M	M	S
500	440	L	M	M	M	M	S
500	450	L	M	M	M	M	S
500	460	M	M	S	M	S	S
500	470	M	M	S	M	S	S

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 75 Degree Frame - Long

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
500	480	M	S	S	M	S	S
500	490	M	S	S	M	S	S
510	440	L	L	L	L	M	M
510	450	L	L	M	L	M	M
510	460	L	L	M	M	M	S
510	470	L	M	M	M	M	S
510	480	L	M	M	M	M	S
510	490	M	M	S	M	M	S
520	450	L	L	L	L	L	M
520	460	L	L	L	L	M	M
520	470	L	L	M	L	M	M
520	480	L	L	M	L	M	M
520	490	L	M	M	L	M	M
530	460	-	L	L	L	L	L
530	470	-	L	L	L	L	M
530	480	-	L	L	L	L	M
530	490	-	L	L	L	L	M
540	470	-	-	L	-	L	L
540	480	-	-	L	-	L	L
540	490	-	-	L	-	L	L

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 85 Degree Frame - Short

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
450	390	M	M	S	M	S	S
450	400	M	M	S	M	S	S
450	410	M	S	S	S	S	S
450	420	S	S	S	S	S	S
460	390	L	M	M	M	M	S
460	400	L	M	M	M	M	S
460	410	M	M	S	M	S	S
460	420	M	M	S	M	S	S
460	430	M	M	S	M	S	S
470	390	L	L	M	L	M	M
470	400	L	L	M	L	M	M
470	410	L	M	M	M	M	S
470	420	L	M	M	M	M	S
470	430	L	M	M	M	M	S
470	440	L	M	M	M	M	S
480	390	L	L	L	L	L	M
480	400	L	L	L	L	L	M
480	410	L	L	L	L	M	M
480	420	L	L	M	L	M	M
480	430	L	L	M	L	M	M
480	440	L	M	M	L	M	M
480	450	L	M	M	L	M	M
490	390	-	L	L	L	L	L
490	400	-	L	L	L	L	L
490	410	-	L	L	L	L	M
490	420	-	L	L	L	L	M
490	430	-	L	L	L	L	M
490	440	-	L	L	L	L	M
490	450	L	L	M	L	L	M
490	460	L	L	M	L	L	M
500	400	-	-	L	-	L	L
500	410	-	-	L	-	L	L
500	420	-	-	L	-	L	L
500	430	-	-	L	-	L	L
500	440	-	-	L	-	L	L
500	450	-	L	L	-	L	L
500	460	-	L	L	L	L	L
500	470	-	L	L	L	L	L
510	410	-	-	-	-	-	L
510	420	-	-	-	-	-	L
510	430	-	-	-	-	-	L
510	440	-	-	-	-	-	L
510	450	-	-	-	-	-	L
510	460	-	-	L	-	L	L

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 85 Degree Frame - Short

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
510	470	-	-	L	-	L	L
510	480	-	-	L	-	L	L
520	420	-	-	-	-	-	-
520	430	-	-	-	-	-	-
520	440	-	-	-	-	-	-
520	450	-	-	-	-	-	-
520	460	-	-	-	-	-	L
520	470	-	-	-	-	-	L
520	480	-	-	-	-	-	L
520	490	-	-	-	-	-	L
530	430	-	-	-	-	-	-
530	440	-	-	-	-	-	-
530	450	-	-	-	-	-	-
530	460	-	-	-	-	-	-
530	470	-	-	-	-	-	-
530	480	-	-	-	-	-	-
530	490	-	-	-	-	-	-
540	440	-	-	-	-	-	-
540	450	-	-	-	-	-	-
540	460	-	-	-	-	-	-
540	470	-	-	-	-	-	-
540	480	-	-	-	-	-	-
540	490	-	-	-	-	-	-

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 85 Degree Frame - Long

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
450	390	S	S	S	S	S	S
450	400	S	S	S	S	S	S
450	410	S	S	S	S	S	S
450	420	S	S	S	S	S	S
460	390	S	S	S	S	S	S
460	400	S	S	S	S	S	S
460	410	S	S	S	S	S	S
460	420	S	S	S	S	S	S
460	430	S	S	S	S	S	S
470	390	M	S	S	S	S	S
470	400	S	S	S	S	S	S
470	410	S	S	S	S	S	S
470	420	S	S	S	S	S	S
470	430	S	S	S	S	S	S
470	440	S	S	S	S	S	S
480	390	M	M	S	S	S	S
480	400	M	M	S	S	S	S
480	410	M	S	S	S	S	S
480	420	M	S	S	S	S	S
480	430	S	S	S	S	S	S
480	440	S	S	S	S	S	S
480	450	S	S	S	S	S	S
490	390	L	M	M	M	S	S
490	400	L	M	M	M	S	S
490	410	L	M	M	M	S	S
490	420	M	M	S	M	S	S
490	430	M	M	S	S	S	S
490	440	M	S	S	S	S	S
490	450	M	S	S	S	S	S
490	460	M	S	S	S	S	S
500	400	L	L	M	M	M	S
500	410	L	L	M	M	M	S
500	420	L	M	M	M	M	S
500	430	L	M	M	M	S	S
500	440	M	M	S	M	S	S
500	450	M	M	S	M	S	S
500	460	M	M	S	M	S	S
500	470	M	S	S	M	S	S
510	410	L	L	L	L	M	M
510	420	L	L	L	L	M	M
510	430	L	L	M	M	M	S
510	440	L	L	M	M	M	S
510	450	L	M	M	M	M	S
510	460	L	M	M	M	M	S

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 85 Degree Frame - Long

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
510	470	M	M	S	M	M	S
510	480	M	M	S	M	M	S
520	420	L	L	L	L	L	M
520	430	L	L	L	L	L	M
520	440	L	L	L	L	M	M
520	450	L	L	M	L	M	M
520	460	L	L	M	L	M	M
520	470	L	L	M	L	M	M
520	480	L	M	M	L	M	M
520	490	L	M	M	M	M	S
530	430	-	L	L	L	L	L
530	440	-	L	L	L	L	M
530	450	-	L	L	L	L	M
530	460	-	L	L	L	L	M
530	470	-	L	L	L	L	M
530	480	L	L	M	L	M	M
530	490	L	L	M	L	M	M
540	440	-	-	L	-	L	L
540	450	-	-	L	-	L	L
540	460	-	-	L	-	L	L
540	470	-	-	L	-	L	L
540	480	-	-	L	L	L	M
540	490	-	L	L	L	L	M