




Rear Wheels - Spoke Guards





1



2



3



4



SCM_269A

Rear Wheels - Spoke Guards

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1564453	Spoke Guard, Transparent (22")	1	2
2	A,B	1564457	Spoke Guard, Transparent (24")	1	2
3	A	1564492	Spoke Guard, Transparent (26")	1	2
4	A,B	1564456	Spoke Guard Kuschall, Black (24")	1	2

NOTE: A - Used with Starec Wheels.
 B - Used with Starec and High Performance Wheels.