





# Rear Wheels - Handrim Hardware



ALUMINUM		BLACK ALUMINUM		TITANIUM		STEEL		SUPERGRIP		SURGE		GEKKO		CURVE	
SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD
1	1	1	1	1	1	1	1	1	1	3	3	1	1	1	1
4	4	4	4	4	4	4	4	4	4	4	6	4	4	4	4
	4		4		4		4		4		7		7		7
KUSCHALL WHEEL															
HIGH PERF. WHEEL															
SPINERGY WHEEL															



ALUMINUM		BLACK ALUMINUM		TITANIUM		STEEL		SUPERGRIP		SURGE	
SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD
3	2			3	2	3	2	3	2		
6	5			6	5	6	5	6	5		
	5				5		5		5		
KUSCHALL WHEEL											
HIGH PERF. WHEEL											
SPINERGY WHEEL											

## Rear Wheels - Handrim Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1495958	Kit, Hardware, Handrim (3 cm)	1	2
2		1495959	Kit, Hardware, Handrim (4 cm)	1	2
3		1495961	Kit, Hardware, Handrim w/ Spokes Guard (4cm)	1	2
4		1496019	Kit, Hardware, Handrim, High Performance (3 cm)	1	2
5		1496020	Kit, Hardware, Handrim, High Performance (4cm)	1	2
6		1497111	Kit, Hardware, Handrim, High Performance w/ Spokes Guard (4cm)	1	2
7	A	1496015	Kit, Hardware, Handrim, Mountain Bike (3cm)	1	2

NOTE: A - GEKKO and CURVE Handrims use this Handrim Hardware.