


Backrest - Adjustable Height Push Handles



SCM_232A

Backrest - Adjustable Height Push Handles

| Item Number | Note | Part Number | Description | Quantity Per | |
|-------------|------|-------------|--|--------------|----------|
| | | | | Package | Assembly |
| 1 | | 1438203 | Kit, Adjustable Height, Integrated Push Handles, Standard, Short | 1 | 1 |
| 2 | | 1549985 | Adjusting Lever, Adjustable Height, Integrated Push Handles (Pair) | 1 | 1 |
| 3 | | 1507734 | Handle, Black (Pair) | 2 | 1 |
| 4 | | 1507752 | Backrest Tube, Straight (265mm), Integrated Push Handles (Pair) | 2 | 1 |
| 5 | | 1507753 | Backrest Tube, Straight (280mm), Integrated Push Handles (Pair) | 2 | 1 |
| 6 | | 1507862 | Backrest Tube, Straight (295mm), Integrated Push Handles (Pair) | 2 | 1 |
| 7 | | 1508438 | Backrest Tube, Straight (310mm), Integrated Push Handles (Pair) | 2 | 1 |
| 8 | | 1508466 | Backrest Tube, Straight (325mm), Integrated Push Handles (Pair) | 2 | 1 |
| 9 | | 1508498 | Backrest Tube, Straight (340mm), Integrated Push Handles (Pair) | 2 | 1 |
| 10 | | 1509103 | Backrest Tube, Straight (355mm), Integrated Push Handles (Pair) | 2 | 1 |
| 11 | | 1509104 | Backrest Tube, Straight (370mm), Integrated Push Handles (Pair) | 2 | 1 |
| 12 | | 1509105 | Backrest Tube, Straight (385mm), Integrated Push Handles (Pair) | 2 | 1 |
| 13 | | 1509106 | Backrest Tube, Straight (400mm), Integrated Push Handles (Pair) | 2 | 1 |
| 14 | | 1509107 | Backrest Tube, Straight (415mm), Integrated Push Handles (Pair) | 2 | 1 |
| 15 | | 1509108 | Backrest Tube, Straight (430mm), Integrated Push Handles (Pair) | 2 | 1 |
| 16 | | 1588217 | Kit, Screw, Backrest Tube | 1 | 1 |