







Chassis - Axle Tube - Aluminum








ANGLE	 CLOTHES GUARD	 SIDEREST	 MUDGUARD	SB34	SB36	SB38	SB40	SB42	SB44	SB46
0°	YES	YES		1	2	3	4	5	6	7
3°	YES	NO		9	10	11	12	13	14	15
	NO	YES		10	11	12	13	14	15	16
6°	YES	NO		18	19	20	21	22	23	24
	NO	YES		19	20	21	22	23	24	25





ANGLE	 CLOTHES GUARD	 SIDEREST	 MUDGUARD		SB34	SB36	SB38	SB40	SB42	SB44	SB46
0°	YES	YES		0°/Vario	1	2	3	4	5	6	7
	YES	NO		3°	1	2	3	4	5	6	7
	NO	YES		3°	2	3	4	5	6	7	8
3°	YES	NO		0°, 3°, Vario	9	10	11	12	13	14	15
	NO	YES		0°, 3°, Vario	10	11	12	13	14	15	16
	YES	NO		6°	10	11	12	13	14	15	16
	NO	YES		6°	11	12	13	14	15	16	17
	YES	NO		3°, 6°	18	19	20	21	22	23	24
6°	NO	YES		3°, 6°	19	20	21	22	23	24	25

Chassis - Axle Tube - Aluminum

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1545827	Axle Tube, Aluminum 0° Complete, 318mm	1	1
2		1545828	Axle Tube, Aluminum 0° Complete, 338mm	1	1
3		1545829	Axle Tube, Aluminum 0° Complete, 358mm	1	1
4		1545830	Axle Tube, Aluminum 0° Complete, 378mm	1	1
5		1545831	Axle Tube, Aluminum 0° Complete, 398mm	1	1
6		1545832	Axle Tube, Aluminum 0° Complete, 418mm	1	1
7		1545833	Axle Tube, Aluminum 0° Complete, 438mm	1	1
8		1545834	Axle Tube, Aluminum 0° Complete, 458mm	1	1
9		1545836	Axle Tube, Aluminum 3° Complete, 324mm	1	1
10		1545837	Axle Tube, Aluminum 3° Complete, 344mm	1	1
11		1545838	Axle Tube, Aluminum 3° Complete, 364mm	1	1
12		1545839	Axle Tube, Aluminum 3° Complete, 384mm	1	1
13		1545840	Axle Tube, Aluminum 3° Complete, 404mm	1	1
14		1545841	Axle Tube, Aluminum 3° Complete, 424mm	1	1
15		1545842	Axle Tube, Aluminum 3° Complete, 444mm	1	1
16		1545843	Axle Tube, Aluminum 3° Complete, 464mm	1	1
17		1545844	Axle Tube, Aluminum 3° Complete, 484mm	1	1
18		1545846	Axle Tube, Aluminum 6° Complete, 354mm	1	1
19		1545847	Axle Tube, Aluminum 6° Complete, 374mm	1	1
20		1545848	Axle Tube, Aluminum 6° Complete, 394mm	1	1
21		1545849	Axle Tube, Aluminum 6° Complete, 414mm	1	1
22		1545850	Axle Tube, Aluminum 6° Complete, 434mm	1	1
23		1545851	Axle Tube, Aluminum 6° Complete, 454mm	1	1
24		1545852	Axle Tube, Aluminum 6° Complete, 474mm	1	1
25		1545853	Axle Tube, Aluminum 6° Complete, 494mm	1	1