


Siderests - Tubular Arm



Siderests - Tubular Arm

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		913.1012	Armrest w/ Standard Cushion, Right	1	1
2		913.1013	Armrest, w/ Standard Cushion, Left	1	1
3		1547610	BoltSet Armrest	1	1
4		1562556	Armrest adapter Right w/ Hardware and Cord	1	1
5		1562557	Armrest adapter Left w/ Hardware and Cord	1	1
6		1562509	Connection Tube 296mm	1	1
7		1562510	Connection Tube 316mm	1	1
8		1562511	Connection Tube 336mm	1	1
9		1562512	Connection Tube 356mm	1	1
10		1562513	Connection Tube 376mm	1	1
11		1562514	Connection Tube 396mm	1	1
12		1562553	Connection Tube 416mm	1	1
13		1554427	Cord, Backrest, Foldable	1	1