

Rear Wheels - Inner Tube and Rim Protection

INNER TUBE

RIM PROTECTOR

22"

RIGHTRUN BLK/GREY

22"

TREADED

24"

RIGHTRUN BLK/GREY

24"

TREADED

24"

MARATHON EVO

24"

SCHWALBE

24"

RIGHTRUN BLUE

24"

MOUNTAIN BIKE

25"

RIGHTRUN BLK/GREY

25"

MARATHON EVO

25"

SCHWALBE

26"

RIGHTRUN BLK/GREY

26"

MARATHON EVO

22"

RIGHTRUN BLK/GREY

22"

TREADED

22"

SOLID TIRE

24"

MOUNTAIN BIKE

24"

RIGHTRUN BLK/GREY

24"

TREADED

24"

MARATHON EVO

24"

SCHWALBE

24"

RIGHTRUN BLUE

24"

SOLID TIRE

24"

RIGHTRUN BLK/GREY

24"

TREADED

24"

MARATHON EVO

24"

SCHWALBE

24"

RIGHTRUN BLUE

24"

SOLID TIRE

25"

RIGHTRUN BLK/GREY

25"

MARATHON EVO

25"

SCHWALBE

26"

RIGHTRUN BLK/GREY

26"

MARATHON EVO

11

Rear Wheels - Inner Tube and Rim Protection

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1570848	Inner Tube (22" x 1")	1	2
2		1570849	Inner Tube (24" x 1")	1	2
3		1570852	Inner Tube, Mountain Bike (24")	1	2
4		1570851	Inner Tube, High-Performance Wheel (26")	1	2
5		1570850	Inner Tube (26" x 1")	1	2
6		1570844	Rim Protection (22")	1	2
7		1570845	Rim Protection (24" Large)	1	2
8		1495962	Rim Protection, High Performance / Spinergy Wheel (24")	1	2
9		1570846	Rim Protection (26")	1	2
11		1571072	Rear Wheel Insert Flat Free (20, 22, 24, 26" Tire)	1	2