


Rear Wheels - Wheel Locks



SCM_339A

Rear Wheels - Wheel Locks

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1562583	Brake, Bent Lever, Left	1	1
2		1562584	Brake, Bent Lever, Right	1	1
3		1562585	Brake, Bent Lever, Long Shoe, Left	1	1
4		1562586	Brake, Bent Lever, Long Shoe, Right	1	1
5		1562587	Brake, Straight Lever, Left	1	1
6		1562588	Brake, Straight Lever, Right	1	1
7		1562589	Brake, Straight Lever, Long Shoe, Left	1	1
8		1562590	Brake, Straight Lever, Long Shoe, Right	1	1
9		1562591	Performance Brake, Left	1	1
10		1562592	Performance Brake, Right	1	1
11		SP1512133	Active Brake Left	1	1
12		SP1512134	Active Brake Right	1	1
13		1562595	Pull Brake, Left	1	1
14		1562596	Pull Brake, Right	1	1
15		1562597	Pull Brake, Long Shoe Left	1	1
16		1562598	Pull Brake, Long Shoe Right	1	1
17		1562599	Brake Clamp, Set	1	1
18		1513177	Brake Handle Standard (Pair)	2	1
19		1513178	Brake Handle, Performance Brake (Pair)	2	1
20		1513179	Brake Lever Extension (Pair)	2	1
21		1579328	Kit, Brake Shoe Standard (40mm)	1	1
22		1579329	Kit, Brake Shoe Long (50mm), 4 Pcs.	1	1