




Options and Accessories - Anti-Tipper



Anti-Tipper Configuration Chart




LEFT SIDE

		SB34	SB36 - SB46
Rear Wheel Size	22 - 25"	7	9
	26"		9



RIGHT SIDE

		SB34	SB36 - SB46
Rear Wheel Size	22 - 25"	8	10
	26"		10



PAIR

		SB34 - SB36		SB38 - SB46	
Rear Wheel Size	22 - 25"	7	8	9	10
	26"			9	10

SCM_341A

Options and Accessories - Anti-Tipper

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1562617	Adaption Holder, Axle Ø30, Complete	1	1
2		1562618	Adaption Holder, Vario-Ax Ø40, Complete	1	1
3		SP1590028	Axle Holder w/ Stop, Adaptation, Ø30	1	1
4		SP1590029	Axle Holder w/ Stop, Adaptation, Ø40	1	1
5		1562619	Kit, Hardware, Adaption Holder	1	1
6		1562621	Kit, Hardware, Adaptions	1	1
7		SP1590076	Anti-Tipper, Short - Left	1	1
8		SP1590077	Anti-Tipper, Short - Right	1	1
9		SP1590078	Anti-Tipper, Long - Left	1	1
10		SP1590079	Anti-Tipper, Long - Right	1	1