


Siderests - Mudguard Sizing



Mudguard Sizing - 22 in. Rear Wheel

Rear Wheel Position (mm)	140	125	110	95	80
SH37	L	L	L	L	M
SH38	L	L	L	M	M
SH39	L	L	M	M	M
SH40	L	M	M	M	M
SH41	M	M	M	M	M
SH42	M	M	M	M	M
SH43	M	M	M	M	M
SH44	M	M	M	M	M
SH45	M	M	M	M	M
SH46	M	M	M	M	M
SH47	M	M	M	M	M
SH48	M	M	M	M	M
SH49	M	M	M	M	M

Siderests - Mudguard Sizing

Mudguard Sizing - 24 in. Rear Wheel

Rear Wheel Position (mm)	140	125	110	95	80
SH37	XL	XL	L	L	L
SH38	XL	L	L	L	L
SH39	L	L	L	L	L
SH40	L	L	L	L	L
SH41	L	L	L	M	L
SH42	L	L	M	M	L
SH43	L	L	M	M	M
SH44	L	M	M	M	M
SH45	M	M	M	M	M
SH46	M	M	M	M	M
SH47	M	M	M	M	M
SH48	M	M	M	M	M
SH49	M	M	M	M	M

Siderests - Mudguard Sizing

Mudguard Sizing - 25 in. Rear Wheel

Rear Wheel Position (mm)	140	125	110	95	80
SH37	XL	XL	XL	XL	L
SH38	XL	XL	XL	L	L
SH39	XL	XL	L	L	L
SH40	XL	L	L	L	L
SH41	L	L	L	M	L
SH42	L	L	M	M	L
SH43	L	L	L	L	M
SH44	L	L	L	M	M
SH45	L	L	M	M	M
SH46	L	M	M	M	M
SH47	M	M	M	M	M
SH48	M	M	M	M	M
SH49	M	M	M	M	M

Siderests - Mudguard Sizing

Mudguard Sizing - 26 in. Rear Wheel

Rear Wheel Position (mm)	140	125	110	95	80
SH37	XL	XL	XL	XL	XL
SH38	XL	XL	XL	XL	XL
SH39	XL	XL	XL	XL	L
SH40	XL	XL	XL	L	L
SH41	XL	XL	L	L	L
SH42	XL	L	L	L	L
SH43	L	L	L	L	L
SH44	L	L	L	L	L
SH45	L	L	L	L	M
SH46	L	L	L	M	M
SH47	L	L	M	M	M
SH48	L	M	M	M	M
SH49	M	M	M	M	M