








Chassis - Axle Tube - Carbon







ANGLE	 CLOTHES GUARD	  SIDEREST MUDGUARD	SB34	SB36	SB38	SB40	SB42	SB44	SB46
0°	YES	YES	1	2	3	4	5	6	7
3°	YES	NO	9	10	11	12	13	14	15
	NO	YES	10	11	12	13	14	15	16
6°	YES	NO	18	19	20	21	22	23	24
	NO	YES	19	20	21	22	23	24	25





ANGLE	 CLOTHES GUARD	  SIDEREST MUDGUARD	 0°/Vario	SB34	SB36	SB38	SB40	SB42	SB44	SB46
0°	YES	YES	0°/Vario	1	2	3	4	5	6	7
	YES	NO	3°	1	2	3	4	5	6	7
	NO	YES	3°	2	3	4	5	6	7	8
3°	YES	NO	0°, 3°, Vario	9	10	11	12	13	14	15
	NO	YES	0°, 3°, Vario	10	11	12	13	14	15	16
	YES	NO	6°	10	11	12	13	14	15	16
	NO	YES	6°	11	12	13	14	15	16	17
	YES	NO	3°, 6°	18	19	20	21	22	23	24
6°	NO	YES	3°, 6°	19	20	21	22	23	24	25

Chassis - Axle Tube - Carbon

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1507565	Axle Tube, Carbon 0° Complete, 318mm	1	1
2		1507566	Axle Tube, Carbon 0° Complete, 338mm	1	1
3		1507567	Axle Tube, Carbon 0° Complete, 358mm	1	1
4		1507568	Axle Tube, Carbon 0° Complete, 378mm	1	1
5		1507569	Axle Tube, Carbon 0° Complete, 398mm	1	1
6		1507570	Axle Tube, Carbon 0° Complete, 418mm	1	1
7		1507571	Axle Tube, Carbon 0° Complete, 438mm	1	1
8		1507572	Axle Tube, Carbon 0° Complete, 458mm	1	1
9		1507575	Axle Tube, Carbon 3° Complete, 324mm	1	1
10		1507576	Axle Tube, Carbon 3° Complete, 344mm	1	1
11		1507577	Axle Tube, Carbon 3° Complete, 364mm	1	1
12		1507578	Axle Tube, Carbon 3° Complete, 384mm	1	1
13		1507579	Axle Tube, Carbon 3° Complete, 404mm	1	1
14		1507580	Axle Tube, Carbon 3° Complete, 424mm	1	1
15		1507581	Axle Tube, Carbon 3° Complete, 444mm	1	1
16		1507582	Axle Tube, Carbon 3° Complete, 464mm	1	1
17		1507583	Axle Tube, Carbon 3° Complete, 484mm	1	1
18		1507585	Axle Tube, Carbon 6° Complete, 354mm	1	1
19		1507586	Axle Tube, Carbon 6° Complete, 374mm	1	1
20		1507587	Axle Tube, Carbon 6° Complete, 394mm	1	1
21		1507588	Axle Tube, Carbon 6° Complete, 414mm	1	1
22		1507589	Axle Tube, Carbon 6° Complete, 434mm	1	1
23		1507590	Axle Tube, Carbon 6° Complete, 454mm	1	1
24		1507591	Axle Tube, Carbon 6° Complete, 474mm	1	1
25		1507592	Axle Tube, Carbon 6° Complete, 494mm	1	1