




Chassis - Axle Tube - Aluminum



0°



3°



				SB36	SB38	SB40	SB42	SB44	SB46	SB48
0°	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	1	2	3	4	5	6	7
	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	2	3	4	5	6	7	8
3°	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	9	10	11	12	13	14	15

Chassis - Axle Tube - Aluminum

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1583478	Axle Tube Aluminum 0° Complete, 348mm	1	1
2		1583479	Axle Tube Aluminum 0° Complete, 368mm	1	1
3		1583480	Axle Tube Aluminum 0° Complete, 388mm	1	1
4		1583481	Axle Tube Aluminum 0° Complete, 408mm	1	1
5		1583482	Axle Tube Aluminum 0° Complete, 428mm	1	1
6		1583483	Axle Tube Aluminum 0° Complete, 448mm	1	1
7		1583484	Axle Tube Aluminum 0° Complete, 468mm	1	1
8		1583637	Axle Tube Aluminum 0° Complete, 488mm	1	1
9		1583486	Axle Tube Aluminum 3° Complete, 368mm	1	1
10		1583487	Axle Tube Aluminum 3° Complete, 388mm	1	1
11		1583488	Axle Tube Aluminum 3° Complete, 408mm	1	1
12		1583489	Axle Tube Aluminum 3° Complete, 428mm	1	1
13		1583490	Axle Tube Aluminum 3° Complete, 448mm	1	1
14		1583491	Axle Tube Aluminum 3° Complete, 468mm	1	1
15		1583492	Axle Tube Aluminum 3° Complete, 488mm	1	1