



Footrest - Adjustable Angle - Adducted Frame



	75°	90°
SHv48	UL40-45	UL38-44
SHv49	UL41-46	UL39-45
SHv50	UL42-47	UL41-46
SHv51	UL43-48	UL42-47

FOOTREST ADJUSTABLE ANGLE ADDUCTED FRAME

SB36	SB38	SB40	SB42	SB44	SB46	SB48
1	1	2	3	3	3	3



4	4	5	6	6	6	6
SB36	SB38	SB40	SB42	SB44	SB46	SB48

Footrest - Adjustable Angle - Adducted Frame

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1583608	Footrest Angle Adjustable 250mm, Complete	1	1
2		1583609	Footrest Angle Adjustable 280mm, Complete	1	1
3		1583610	Footrest Angle Adjustable 308mm, Complete	1	1
4		1583613	Footplate, 250mm w/ Clamping Part	1	1
5		1583614	Footplate, 280mm w/ Clamping Part	1	1
6		1583615	Footplate, 308mm w/ Clamping Part	1	1
7		1583616	Kit, Telescopic Footrest, Adjustable Angle (Pair)	2	1
8		1583602	Footrest Sleeve, Aluminum	1	1