










# Electronics & Pendants - BP1180/76



SICC\_123A



## Electronics & Pendants - BP1180/76

| Item Number                                                                                                                                                                                                         | Note | Part Number     | Description                                                                                  | Quantity Per |          |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------|-----------------|----------------------------------------------------------------------------------------------|--------------|----------|
|                                                                                                                                                                                                                     |      |                 |                                                                                              | Package      | Assembly |
| 1                                                                                                                                                                                                                   |      | MTRM2133-NLA    | NLA - Actuator, Head (Red Plug) rEF.<br>#61832-1005270-01                                    | 1            | 1        |
| 2                                                                                                                                                                                                                   |      | KIT80-NLA       | NLA - Kit, Bracket Trunion for Head Motor                                                    | 1            | 1        |
|                                                                                                                                                                                                                     |      |                 | Pin, Cotter Bow-Tie (1/4 x 3/32 x 5/16")                                                     | 1            | 1        |
|                                                                                                                                                                                                                     |      |                 | Screw, Torx HD Taptite LG (8-32 x 1-1/4")                                                    | 1            | 4        |
|                                                                                                                                                                                                                     |      |                 | Bracket, Trunnion Plated                                                                     | 1            | 1        |
|                                                                                                                                                                                                                     |      |                 | Clevis Pin (3/8 x 2-1/2" L)                                                                  | 1            | 2        |
|                                                                                                                                                                                                                     |      |                 | NLA - Wrap, Tie (8")                                                                         | 1            | 8        |
|                                                                                                                                                                                                                     |      |                 | Sheet, Instruction, Bracket Solo Head/HiLo Motors                                            | 1            | 1        |
| 3                                                                                                                                                                                                                   | A    | MTRM2134-NLA    | NLA - Foot Actuator (Blue Plug) Ref.<br>#61831-0625270-01                                    | 1            | 1        |
| 4                                                                                                                                                                                                                   |      | MTRM2378        | HiLo Actuator (White Plug) Ref. #AP-8390-0114                                                | 1            | 1        |
| 5                                                                                                                                                                                                                   |      | JNCTNBOX-NLA    | NLA - Junction Control Box                                                                   | 1            | 1        |
| 6                                                                                                                                                                                                                   |      | JNCTNBOXMEM-NLA | NLA - Junction Box Membrane                                                                  | 1            | 1        |
| 7                                                                                                                                                                                                                   | B    | 1190227         | Pendant, Six Function (w/ Loop Adhesive Back Strip and DB9 Connector) Linak # PNDNT-6FNC-1-2 | 1            | 1        |
| 8                                                                                                                                                                                                                   | B    | 1190232         | Pendant, Six Function (w/ Strap and DB9 Connector) Linak # PNDNT-6FNC-1-2                    | 1            | 1        |
| 9                                                                                                                                                                                                                   | C    | PNDNT6FNCIV3    | Pendant, Six Function w/ RJ48 Connector) Linak # PNDNT-6FNC-1-2                              | 1            | 1        |
| 10                                                                                                                                                                                                                  |      | 1190225         | Cable, Conversion                                                                            | 1            | 1        |
| 11                                                                                                                                                                                                                  |      | 1186847         | Holster, Pendant                                                                             | 1            | 1        |
| 12                                                                                                                                                                                                                  | B    | PNDNTYCABLE     | Cable, Y, Pendant                                                                            | 1            | 1        |
| 13                                                                                                                                                                                                                  | C    | JNCTNDC40       | Cable, Y, Pendant                                                                            | 1            | 1        |
| 14                                                                                                                                                                                                                  | C    | PNDNTEXTCABLE   | Cable, Pendant                                                                               | 1            | 1        |
| 15                                                                                                                                                                                                                  |      | JNCTNDC20C      | Membrane Pad for 151 Only                                                                    | 1            | 1        |
| NOTE: A - Foot Actuator Kit (KIT84) is No Longer Available for Beds built prior to 12/15/08. Part #MTRM2134 used for Brds built as of 12/08.<br>B - Beds built prior to 11/22/07.<br>C - Beds built after 11/22/07. |      |                 |                                                                                              |              |          |