


Steel and Frame Parts - BS3180/76



Steel and Frame Parts - BS3180/76

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1185149-NLA	NLA - Head Gatch Assembly (Black Insert w/ Steel Pin)	1	1
2		GATFOOTASSY 1-NLA	NLA - Foot Gatch (Black Insert with Steel Pin)	1	1
3	A	HILOUPPERAS SY-NLA	NLA - Upper Gatch "4Arm" (Black Insert w/ Steel Pin)	1	1
4	A	HILOLOWERAS Y1-NLA	NLA - Gatch, Lower 8 Arm Head	1	1
5	A	HILOLOWERAS Y2-NLA	NLA - Gatch, Lower 6 Arm Foot	1	1
6		DRVTUBESOLO P1-NLA	NLA - HiLo Drive Tube	1	1
7		LNK0017PL-NLA	NLA - Linkage, Drive Tube	1	1
8		TBR78161812P N-NLA	NLA - Arm, Solo Control	1	1
9	A	LNK0013PNT-N LA	NLA - Linkage, Hilo Head/Foot	1	1
10		LNK0018PL-NLA	NLA - HiLo HD/FT Leg Linkage	1	1
11	A	WHLBRACKASS Y3-NLA	NLA - Caster Assembly (Non-Lock)	1	1
12	A	WHLBRACKASS Y4-NLA	NLA - Caster Assembly, Lock	1	1
13		MNT0018GPNT- NLA	NLA - Crank Mount Plate	1	1
14		LWRFRASSY2- NLA	NLA - Lower Frame	1	1
15		LEG0003DASSY 1	Leg, "H" (Head/Foot)	1	1
16		LEG0004DASSY -NLA	NLA - Head End Caster Leg Assembly	1	1
17		LEG0003DASSY -NLA	NLA - Foot End Caster Leg Assembly	1	1
18		1186525	Post, Head	1	1
19		1186524	Post, Foot	1	1
20		KIT78	Kit, Head & Foot Board Mounting Hardware	1	1
			Knob, Three Arm	1	2
21		GLDBEDFOOT	Foot, Bed - Rubber Santoprene	1	1
22		1186862	Plug, Under Bed Foot	1	1
23		PIN38412-NLA	NLA - Pin, 4-1/2" Straight	1	1
24		PIN380212	Clevis Pin (3/8 x 2-1/2" L)	1	1
25		PINBOW516-NL A	NLA - Pin, Cotter, Bowtie	1	1
26		MATTKPRKIT	Kit, Mattress Keeper, Left and Right Corners	1	1
			Mattress Keeper, Left Corner	1	2
			Mattress Keeper, Right Corner	1	2
27	B	IHECTA-NLA	NLA - Trapeze Adapter	1	1
28		FRENCHERMAL E-NLA	NLA - Plate, Frencher - Male	1	1
29	D	GLDBUMPER	Rubber Bed Foot & Wall Bumper (Optional, 2 per Wireform)	1	1

Steel and Frame Parts - BS3180/76

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
30	A,C	KIT19EARRO-NLA	NLA - Extension Kit - 4"	1	1
31		WBBSOLO76P-NLA	NLA - Wireform, 76"	1	1
32		WBBSOLO80P-NLA	NLA - Wireform, 80"	1	1
NOTE: A - Order through SFI Medical Solutions B - Trapeze Adapter comes with Frame, Mounting Hardware and (2 ea. Male Frencher Plates) C - Does not come with Bolster Pad, must use 84" Mattress. D - Optional (2 per Wireform)					