



Electronics & Pendants


NOTE: This Bed is obsolete and only parts listed are available to order.
NOTE: All Frame Parts are obsolete and No Longer Available. For Deck Replacements, see Echo 2380/76 Deck Replacement Guide.



Shaft Approx. 3 Ft. - Red Plug





Shaft Approx. 1 Ft. - Blue Plug




Shaft Approx. 1-1/2 Ft. - White Plug



Comes w/ Groundwire and Clips



DB9 Connector



DB9 Connector

SICC_110A

Electronics & Pendants

NOTE: This Bed is obsolete and only parts listed are available to order.

NOTE: All Frame Parts are obsolete and No Longer Available. For Deck Replacements, see Echo 2380/76 Deck Replacement Guide.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	MTRM1703	Head Actuator (Red Plug)	1	1
2	B	MTRM1702	Foot Actuator (Blue Plug)	1	1
3	C	MTRM1747	NLA - Hilo Actuator (White Plug)	1	1
4	D	JNCTNBOXKIT1 -NLA	NLA - Junction Control Box (Comes w/ Ground Wire & Clips)	1	1
5	E	1190233	Kit, Pendant, Full/Partial Rails, Linak # PNDNT-6FNC-1-2 (With DB9 Connector)	1	1
6	E	PNDNTKIT3-NL A	NLA - Kit, Pendant, Assist Rails (Old Square), Linak # PNDNT-6FNC-1-2 (With DB9 Connector) (May 01 - July 02)	1	1

NOTE: A - Will work on old Allto Beds. VO7508AW76U
 B - Will work on old Allto Beds. VO7508AX76U
 C - Reference 6183M-1905220-01
 D - J Box Prior to 11/29/01. NAX5105-3A29MPZ
 E - These Pendants are not designed for the Allto, but these will work. However the Pendant cannot be fastened to the Control Box using the screws.