



Caster Kits - IH6065 Series

Models 6065A, 6065DA, 6065PT & 6065WD



SPT_027A

Caster Kits - IH6065 Series

Models 6065A, 6065DA, 6065PT & 6065WD

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1126605	Kit, Casters - Front	1	1
2	B	1126606	Kit, Casters - Rear	1	1
NOTE: A - Kit includes: (2) Casters, (2) Pan Head Stem Bolts (3/8-16 x 2-3/8) and (2) Nylock Nuts (3/8-16). B - Kit includes: (2) Locking Casters, (2) Hex Head Stem Bolts (1/4-20 x 1-1/2), (4) Washers (6.5 mm) and (2) Hex Nuts (1/4-20).					