



Caster Kits - IH6077 Series

Models 6077A

①



②



Caster Kits - IH6077 Series Models 6077A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1192445	Kit, Casters - Front	1	1
2	B	1126606	Kit, Casters - Rear	1	1
NOTE: A - Kit includes: (2) Casters, (2) Pan Head Stem Bolts (1/4-20 x 1-1/2), (2) Nylock Nuts (1/4-20), and (4) Washers (6.5mm). B - Kit includes: (2) Locking Casters, (2) Hex Head Stem Bolts (1/4-20 x 1-1/2), (4) Washers (6.5mm) and (2) Hex Nuts (1/4-20).					