



Tray and Side Panel Kits - IH6074 Series Models 6074A



SPT_031A

Tray and Side Panel Kits - IH6074 Series Models 6074A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1126692	Tray Assembly	1	1
2	B	1126611-NLA	NLA - Kit, Side Panels - Plastic (31.8" x 16.8")	1	1
			Panel, Side - Vacuum Formed - Right	1	1
			Panel, Side - Vacuum Formed - Left	1	1
2	B	1164310	Kit, Side Panels - Plastic (32.33" x 17.24")	1	1
			Panel, Side - Vacuum Formed - Right	1	1
			Panel, Side - Vacuum Formed - Left	1	1
3		1134780	Kit, Tray Maintenance	1	1

NOTE: A - Includes: Tray Assembly with (2) Mounting Tubes and (2) Phillips Head Screws (M6 x 34).
 B - To determine which part number to use, measure the outside dimension (length and width) to match part number to order.