

Tray and Side Panel Kits - IH6077 Series Models 6077A



SPT_032A

Tray and Side Panel Kits - IH6077 Series Models 6077A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	C,A,D	1126728-NLA	NLA - Kit, Side Panels - Plastic (32-1/2" x 18-1/4")	1	1
1	E,F	1165090-NLA	NLA - Kit, Side Panels, Right/Left (IH6077A)	1	1
			Panel, Side, Vacuum Formed - Right	1	1
			Panel, Side, Vacuum Formed - Left	1	1
2	B	1126727	Tray Assemblies with Mounting Tubes	1	1
NOTE: A - Kit includes: (2) Side Panels. B - Kit includes: (2) Tray Assemblies with Mounting Tubes. C - To determine which part number to use, measure the outside dimension (length and width) to match part number to order. D - Used on Recliner Models built before October 2008. E - Used on Recliner Models built after October 2008. F - Kit includes Left and Right Side Panels.					