



# Upholstery Kits - IH6065 Series


Models 6065A, 6065DA, 6065PT & 6065WD



①



②



③

SPT\_032A

## Upholstery Kits - IH6065 Series

### Models 6065A, 6065DA, 6065PT & 6065WD

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	Back and Headrest Upholstery Kits	1	1
2	B	See Chart	Arm Rest Upholstery Kits	1	1
3	C	See Chart	Seats with Footrest Pad and Skirt Upholstery Kits	1	1
NOTE: A - Kit includes: Back, Headrest, (8) Hex Head Screws (1/4-20 x 1), (8) Washers (6.5mm) and (2) Angle Brackets. B - Kit includes: (2) Arm Rests and (8) Phillips Head Screws (1/4-20 x 5/8). C - Kit includes: Seat with Foot Rest, Skirt and (8) Phillips Head Screws (1/4-20 x 1).					

## Upholstery Kits - IH6065 Series

### Models 6065A, 6065DA, 6065PT & 6065WD

**Chart 1 - Back and Headrest**

Color	IH6065A	IH6065DA	IH6065PT	IH6065WD
IH60 - Rosewood	1126589	1126612-NLA	1126636-NLA	1126652-NLA
IH61 - Blue Ridge	1126590	1126613-NLA	1126637-NLA	1126653
IH68 - Jade	1126591-NLA	1126614-NLA	1126638-NLA	1126654-NLA

**Chart 2 - Arm Rest**

Color	IH6065A, IH6065DA, IH6065PT, IH6065WD
IH60 - Rosewood	1126581
IH61 - Blue Ridge	1126582
IH68 - Jade	1126583

**Chart 3 - Seat with Footrest Pad and Skirt**

Color	IH6065A	IH6065DA	IH6065PT	IH6065WD
IH60 - Rosewood	1126597	1126620-NLA	1126644-NLA	1126660
IH61 - Blue Ridge	1126598	1126621-NLA	1126645-NLA	1126661
IH68 - Jade	1126599	1126622-NLA	1126646-NLA	1126662