



Upholstery Kits - IH6074 Series

Models 6074A



SPT_034A

Upholstery Kits - IH6074 Series Models 6074A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	Back and Headrest Upholstery Kits	1	1
2	B	See Chart	Arm Rest Upholstery Kits	1	1
3	C	See Chart	Seat with Footrest Upholstery Kits	1	1
NOTE: A - Kit includes: Back, Headrest, (8) Hex Head Screws (1/4-20 x 1), (8) Washers, (2) Hex Nuts (1/4-20) and (2) Angle Brackets. B - Kit includes: (2) Arm Rests and (8) Phillips Head Screws (1/4-20 x 5/8). C - Kit includes: Seat with Foot Rest, Skirt and (8) Phillips Head Screws (1/4-20 X 1).					

Upholstery Kits - IH6074 Series

Models 6074A

Chart 1 - Back and Headrest Upholstery Kits

Color	Kit Number
IH60 - Rosewood	1126676
IH61 - Blue Ridge	1126677
IH68 - Jade	1126678

Chart 2 - Arm Rest Upholstery Kits

Color	Kit Number
IH60 - Rosewood	1126668
IH61 - Blue Ridge	1126669
IH68 - Jade	1126670

Chart 3 - Seat with Footrests Upholstery Kits

Color	Kit Number
IH60 - Rosewood	1126684
IH61 - Blue Ridge	1126685
IH68 - Jade	1126686-NLA