



Upholstery Kits - IH6077 Series Models 6077A



SPT_035A

Upholstery Kits - IH6077 Series Models 6077A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	Back and Headrest Upholstery Kits	1	1
2	B	See Chart	Arm Rest Upholstery Kits	1	1
3	C	See Chart	Seat with Footrest Upholstery Kits	1	1
4	D	See Chart	Wing Pad Kits	1	1
<p>NOTE: A - Kit includes: Back, Headrest, (8) Hex Head Screws (1/4-20 X 1), (2) Hex Nuts (1/4-20), (8) Washers (6.5mm), and (2) Angle Brackets. B - Kit includes: (2) Arm Rests and (8) Phillips Head Screws (1/4-20 x 5/8). C - Kit includes: Seat with Foot Rest, Skirt and (4) Phillips Head Screws (1/4-20 x 1). D - Kit includes: (2) Wing Pads with Backs and (4) Phillips Head Screws (1/4-20 x 1).</p>					

Upholstery Kits - IH6077 Series

Models 6077A

Chart 1 - Back and Headrest Upholstery Kits

Upholstery Color	Kit Number
IH60 - Rosewood	1126701
IH61 - Blue Ridge	1126702
IH68 - Jade	1126703-NLA

Chart 2 - Arm Rest Upholstery Kits

Upholstery Color	Kit Number
IH60 - Rosewood	1126693
IH61 - Blue Ridge	1126694
IH68 - Jade	1126695

Upholstery Kits - IH6077 Series

Models 6077A

Chart 3 - Seat with Footrests Upholstery Kits

Upholstery Color	Kit Number
IH60 - Rosewood	1126709-NLA
IH61 - Blue Ridge	1126710
IH68 - Jade	1126711-NLA

Chart 4 - Wing Pads Upholstery Kits

Upholstery Color	Kit Number
IH60 - Rosewood	1126717
IH61 - Blue Ridge	1126718
IH68 - Jade	1126719