




Electronics and Pendants - BC2380/76



Red Plug



Shaft Approx. 1 Ft. - Blue Plug




Shaft Approx. 1-1/2 Ft. - White Plug




Comes w/ Ground Wire and Clips




DB9 Connector




DB9 Connector



DB9 Connector



DB9 Connector



SICC_098A

Electronics and Pendants - BC2380/76

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		MTRM1633	Actuator, Head (Red Plug) Ref. #61831-0795280-01	1	1
1		MTRM1633HC8 0-NLA	NLA - Actuator, Head - Bed (80")	1	1
2		MISC022	Cable, Extension - Head Motor	1	1
3		MTRM1702	Foot Actuator (Blue Plug)	1	1
4		MTRM1750-NLA	NLA - Hilo Actuator (White Plug)	1	1
5	A	JNCTNBOX-NLA	NLA - Junction Control Box	1	1
6	B	JNCTNBOXKIT1 -NLA	NLA - Junction Control Box (Comes w/ Ground Wire & Clips)	1	1
7		PNDNT2FNCIV	Pendant, 2 Function	1	1
8	C	1190228	Kit, Pendant, New Assist Rails, Linak # PNDNT-6FNC-1-2	1	1
9	C	1190233	Kit, Pendant, Full/Partial Rails, Linak # PNDNT-6FNC-1-2 (With DB9 Connector)	1	1
10	C	PNDNTKIT3-NL A	NLA - Kit, Pendant, Assist Rails (Old Square), Linak # PNDNT-6FNC-1-2 (With DB9 Connector) (May 01 - July 02)	1	1
11		PNDNTYCABLE	Cable, Y, Pendant	1	1
12		1190225	Cable, Conversion	1	1
NOTE: A - Beds manufactured After 11/29/01 B - Beds manufactured Before 11/29/01 C - Reference Linak #PNDNT-6FNC-1-2					