





Electronics and Pendants - BC1180/76




Red Plug




Shaft Approx. 1 Ft. - Blue Plug





Shaft Approx. 4-1/2 Ft. - White Plug



DB9 Connector




DB9 Connector




DB9 Connector




9




RJ48 Connector




11



12



13



14

Electronics and Pendants - BC1180/76

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		MTRM2133-NLA	NLA - Actuator, Head (Red Plug) rEF. #61832-1005270-01	1	1
2		MISC022	Cable, Extension - Head Motor	1	1
3		MTRM2102-NLA	NLA - Actuator, Foot (Blue Plug) 115V ~ 60Hz	1	1
4		MTRM2150	Actuator, HILO (White Plug) 115V ~ 60Hz	1	1
5		JNCTNBOX-NLA	NLA - Junction Control Box	1	1
6	A	PNDNT2FNCIV	Pendant, 2 Function	1	1
7	A	1190227	Pendant, Six Function (w/ Loop Adhesive Back Strip and DB9 Connector) Linak # PNDNT-6FNC-1-2	1	1
8	A	1190232	Pendant, Six Function (w/ Strap and DB9 Connector) Linak # PNDNT-6FNC-1-2	1	1
9	A	1190225	Cable, Conversion	1	1
10	B	PNDNT6FNCIV3	Pendant, Six Function w/ RJ48 Connector) Linak # PNDNT-6FNC-1-2	1	1
11		1186847	Holster, Pendant	1	1
12	A	PNDNTYCABLE	Cable, Y, Pendant	1	1
13	B	JNCTNDC40	Cable, Y, Pendant	1	1
14	B	PNDNTEXTCABLE	Cable, Pendant	1	1

NOTE: A - Beds manufactured Before 11/22/07
B - Beds manufactured After 11/21/07